MTC306

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


EN

MTC306

Additional information is available on the manufacturer's website.

Introduction

The probe is a combined silver and potassium nitrate electrode with refillable electrolyte, a glass body and an integrated temperature sensor. The probe is used for halides determination in potentiometric titrations.

Warning

⚠ Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Items to collect:

AT1000 series titrator, deionized water, no-lint cloth, filling solution

1 Connect the probe

2 Calibration

⚠ Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.

3 Measurement

Technical information

Range	±1200 mV
Temperature	0 to 80 °C (32 to 176 °F)
Filling solution	KNO ₃ 1M
Storage solution	KNO ₃ 1M

Maintenance and storage

Soak the probe in a cleaning solution. Rinse with deionized water. If necessary, rub the silver rod with polishing bands to remove contamination. Soak in acetone and rinse with distilled water.

Add filling solution to the probe when the filling solution level is low.

For storage, keep the electrode in the soaker bottle with storage solution.

Warranty

1 year for manufacturing faults. This warranty does not cover improper use or wear.

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