DOC022.98.00743

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Section 1

Description of the piston pipette


Figure 1 Description of the piston pipette

1	Piston pipette 0.1 to 1 mL	5	Tip ejector
2	Push button	6	Cone
3	Ejector button	7	Piston pipette 1 to 5 mL
4	Volume scale		

Description of the piston pipette

The piston pipette is a handy, general-purpose microliter pipette for transferring exact volumes of liquid. It operates on the principle of air expulsion, making use of disposable tips. The volume that is to be pipetted is shown in a clearly visible display on the pipette handle. All piston pipettes are fitted with a mechanical tip ejector system.

1.1 Pipette tips

We recommend the use of original HACH pipette tips when working with the piston pipette. The pipette tips are made of polypropylene and can be autoclaved at 121 °C/250 °F.

Section 2 Using the pipette

2.1 Setting the volume

The pipette volume is set with the push button at the end of the handle. Turn the button anticlockwise to increase the volume and clockwise to decrease it. Make sure that the volume setting is correct (clicks into place) and the volume shown in the display is completely visible.

Figure 2 Setting the volume

Important Note: *The given precision and accuracy only apply for the volume range specified on the pipette. Do not set the volume to a value outside of this range. Do not use force to turn the push button to a volume that is outside of the specified range. This could damage the pipette.*

2.2 Tip ejection

Each piston pipette is fitted with a tip ejection system to prevent possible contamination. Press thumb down on ejector button and remove tip from tip cone.

Using the pipette

2.3 Pipetting techniques

The pipette is operated by means of the push button at the end of the handle.

Maximum accuracy can be achieved by carrying out the following steps:

1. Move the piston up or down slowly and evenly (especially when pipetting highly viscous liquids).
2. Do not allow the push button to spring back into its original position.
3. Use light pressure to fit a new tip to the cone, turning the tip slightly as you do so.
4. Make sure that the tip is firmly in position on the tip cone and there are no foreign bodies between the tip and the cone.
5. Wet the tips with the liquid by repeatedly filling and emptying them.
6. Hold the pipette almost vertical (no more than 10° from the vertical) while liquid is being drawn into the tip.
7. The temperature of the pipette and the tip should match that of the liquid.

2.4 Standard technique

1. Press the push button to the first pressure point.
2. Dip the pipette tip under the surface of the liquid (2–3 mm) and slowly release the push button. Take the tip out of the liquid and touch it against the inside wall of the vessel to remove excess liquid.
3. Place the tip against the inside wall of the vessel and press lightly and evenly on the push button until the first pressure point is reached. Hold the pipette in this position. After about one second, press the push button until the

second pressure point is reached. This completely empties the tip.

4. Allow the push button to return to its rest position.


Figure 3 Standard technique

a Rest position	c Second pressure point
b First pressure point	

Important Note: After use change pipette tip!

Section 3 Maintenance

3.1 Daily check

Check the pipette daily for dust and external soiling before starting work or after finishing work. Pay special attention to the tip cone. Clean the pipette with 70% alcohol. Do not use any other solvents to clean the pipette.

3.2 Maintenance

If the pipette is used each day, maintenance should be carried out at least **twice each year**. Contact the Hach Technical support.

Section 4 Calibration

The pipette is calibrated in the factory with distilled water at 22 °C/72 °F. No recalibration is needed for routine use. The pipette can, however, be recalibrated for other temperatures and solutions of different viscosities. An appropriately sensitive analytical balance, a small beaker and distilled water are needed to carry out the calibration.

4.1 Checking the calibration

1. Set the pipette volume as follows in [Table 1 Checking the calibration](#):

Table 1 Checking the calibration

Pipette (μL)	Volume setting	permissible weight range (mg)
100	100	99.2–100.8
200	200	197.5–202.5
500	500	497.0–503.0
1000	1000	994–1006
2000	2000	1990–2010
5000	5000	4975–5025
100–1000	300	298–302
1000–5000	2000	1990–2010

2. Use light pressure to place a tip on the tip cone, turning the tip slightly.
3. Pipette distilled water into previously weighed beaker at least **five times**. Determine each weight in mg to the first decimal place and compare this with the weight shown above. If the weight is not within this range, recalibrate the pipette.

4.2 Recalibration

The calibration is carried out with the supplied calibration key.

1. Place the calibration key in the holes of the calibration groove at the top end of the handle.
2. Increase the volume by turning in clockwise direction or reduce it by turning in anticlockwise direction.

Figure 4 Recalibration

1	Calibration key
---	-----------------

3. Afterwards check the volume as described above
[Section 4 Calibration on page 11](#).

Important Note: The pipette can only be calibrated within the limits of the volume range for which it is meant to be used.

Section 5 Troubleshooting

The following table is intended to help you to identify and eliminate errors:

Table 2 Troubleshooting

Error	Possible cause	Response
Leakage	Tip not attached properly	Fix tip firmly in place
	Foreign bodies between cone and tip	Clean the cone and attach a new tip
Inaccurate pipetting	Incorrect operation	Use pipette as described in the manual
	Tip not firmly attached	Attach tip firmly
	Calibration has changed, e.g. due to faulty operation	Recalibrate

Section 6 **Warranty, liability and complaints**

The manufacturer warrants that the product supplied is free of material and manufacturing defects and undertakes the obligation to repair or replace any defective parts at zero cost.

The warranty period for these pipettes is 12 months (EU 24 months).

With the exclusion of the further claims, the supplier is liable for defects including the lack of assured properties as follows: all those parts that can be demonstrated to have become unusable or that can only be used with significant limitations due to a situation present prior to the transfer of risk, in particular due to incorrect design, poor materials or inadequate finish will be improved or replaced, at the supplier's discretion. The identification of such defects must be notified to the supplier in writing without delay, however at the latest 7 days after the identification of the fault. If the customer fails to notify the supplier, the product is considered approved despite the defect. Further liability for any direct or indirect damages is not accepted.

If specific maintenance work defined by the supplier is to be performed within the warranty period by the customer (maintenance) and these requirements are not met, claims for damages due to the failure to comply with the requirements are rendered void.

Any further claims, in particular claims for consequential damages cannot be made.

Consumables and damage caused by improper handling, poor installation or incorrect use are excluded from this clause.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)