
LANGE 

DOC022.98.80410

QuikChem[®] 8500 Series 2

09/2013, Edition 1

Maintenance and Troubleshooting
Wartung und Fehlerbehebung
Maintenance et dépannage
Mantenimiento y solución de problemas

English 3

Deutsch 20

Français 39

Español 58

Table of contents

[Additional information](#) on page 3

[General information](#) on page 3

[Maintenance](#) on page 4

[Troubleshooting](#) on page 14

Additional information

Additional information is available on the manufacturer's website.

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. <i>Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.</i>
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates the need for protective eye wear.
	This symbol indicates that the marked item can be hot and should not be touched without care.
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.
	This symbol indicates that the object is heavy.

Maintenance

⚠ DANGER	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.
⚠ DANGER	
	Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.
⚠ WARNING	
 	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Maintenance schedule

Table 1 shows the recommended schedule of maintenance tasks. Facility requirements and operating conditions can increase the frequency of some tasks.

Table 1 Maintenance schedule

Component	Procedure	As necessary	Daily	Weekly	Monthly	6 months	Yearly
Sampler	Clean surfaces	X	X				
	Clean rods/moving parts				X		
	Remove loose particles off the y-axis lead screw with a dry, lint-free cloth.			X			
	Replace the pump tubes on the rinse station assembly.				X		
Dilutor	Clean surfaces	X	X				
	Prime dilutor with DI water before another diluent is used		X				
	Replace all tubing						X
Pump	Spray silicone on cloth and rub on rollers				X		
	Replace pump tubes				X		
	Clean surfaces	X	X				
	Rinse cartridges		X				
	Clean pump tube adapters				X		
Injection valves	Clean ports and valve connectors			X			
Manifolds	Clean unions and tees				X		
	Replace O-rings					X	
	Replace all tubing						X
Detectors	Dry and clean all surfaces	X	X				
Flow cells	Replace flares and O-rings						X
System unit	Keep dry and clean		X				
	Replace the heater mandrel	X					
Computer	Clean hard drive						X

How to clean spills

▲ CAUTION	
	Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

1. Obey all facility safety protocols for spill control.
2. Discard the waste according to applicable regulations.

Clean the sampler

Clean off spills immediately. Neutralize acid spills with a commercial acid neutralizer solution.

Clean the sample tray and sampler core unit with a lint-free cloth dampened with a lab-grade cleaning agent. Do not let the cleaning agent come in contact with the lead screws. Remove the cleaning agent with a lint-free towel dampened with water.

Clean the modules

⚠ DANGER	
	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Clean the exterior of the modules with water and a paper towel. Remove chemical spills according to MSDS instructions.

Shut down the system

1. To use a manifold with a column, such as the cadmium column for nitrate, make the column offline. Set the switching valve to the bypass position. If the manifold has a column, but does not have a switching valve, set the pump power to off and replace the column with a piece of Teflon tubing. Start the pump.

Note: Some methods make it necessary for some reagents (typically, the buffer) to be removed last. Look at the QuikChem method for instructions.

2. Remove the reagent tubing from each reagent and rinse the tubing and glass weights. Keep the lines and weights in rinse solutions.
3. If the QuikChem method recommends a rinse solution, install all reagent transmission tubing in the solution and operate the pump for five minutes at standard speed.

Note: This is a critical step in the preventive maintenance of the manifolds.

4. Install the transmission tubing in deionized water and let the system rinse for 10 minutes at standard speed.
5. If the manifold will not be used for more than three days or it is removed and kept in storage:
 - a. Remove the transmission tubing from the deionized water.
 - b. Use the pump to remove all liquid from the manifold.
6. Set the power to the pump to off. Release the pump tube cartridges.
7. Close all files on the computer.
8. Set the power to the power strip to off.

Note: If a lower set-point temperature is not given in the software, the heater controller stays at the original set temperature. If the method makes it necessary for a temperature higher than 60 °C (140 °F), lower the temperature setting to 60 °C (140 °F) or less before power is removed from the system.

Heater mandrel maintenance

⚠ WARNING	
	Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

⚠ CAUTION


Burn hazard. Turn off power and allow the instrument to cool before this procedure.

To replace a mandrel:

1. Cut the length of tubing specified in the QuikChem method. Refer to [Figure 1](#).
2. Replace the tubing and tape on the mandrel. Refer to [Figure 2](#).


Figure 1 Tubing preparation

1 Tubing	4 Tubing for other connections ¹
2 Thermal tape, high temperature	5 Start and stop marks
3 Thermal tape, heater	6 Tubing winds around the lamp/heater block

¹ For lengths and other information, refer to the QuikChem method.

Figure 2 Wind a heater mandrel

Replace the dilutor tubing

Replace all tubes and the seals of the barrel cavity once a year. Contact technical support for assistance.

When a diluent other than deionized water (e.g., KCl solution) is used, prime the dilutor with deionized water at least two times, immediately after use. Make sure that the diluent is flushed out completely.

1. In the Omnion software, go to Configuration>Autosamplers.
2. Push the **PRIME DILUTOR** button.
3. If the dilutor leaks from the connectors, replace the Teflon tubing. If the barrel leaks, contact technical support for assistance.
Note: Always degas the diluent with helium.
4. Keep all surfaces clean. Clean the exterior of the modules with water and a paper towel. Dry all surfaces thoroughly. If an acid spill occurs, neutralize the acid with baking soda and then clean all the surfaces.

Sampler maintenance

⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

NOTICE

Never lubricate the leads screws. The lead screw nuts are compounded with a dry film lubricant. If oil is used on the lead screws, it will cause gumming, galling and binding of the sample probe assembly.

Clean and examine the pump

The pump is susceptible to acid damage. The cartridges and cartridge holders can be damaged due to acid. Do the steps that follow to clean the pump.

1. After daily use, rinse the cartridges in deionized water to wash off spills that come in contact with them.
2. Clean the pump surfaces, except the rollers, with a wet cloth. Dry all surfaces well.
3. Remove rust from the rollers and clean with steel wool. Apply a very light coat of silicone spray to a lint-free cloth and hold it on the moving rollers. The silicone will act as a rust inhibitor.
4. Look for wear, cracks or acid damage on the cartridges and holders. Replace as necessary. Make sure that the flow is correct at the time of analysis.
5. Replace pump tubing as soon as it starts to show signs of wear, such as sections that are pushed flat (approximately once per month). If the usage is high, then replace the pump more often.
6. If a pump tube bursts, immediately clean all cartridges, holders and the pump. If a spill is not cleaned immediately, it can cause damage to the pump and/or parts.
7. If the pump speed looks too fast or too slow when the speed is set at 35, do the pump test to make sure that the pump operates correctly. Refer to the expanded user manual available on the manufacturer's website.
8. Keep all surfaces clean. Use a damp cloth to clean the module surfaces. Dry all surfaces thoroughly. If an acid spill occurs, use baking soda to make the acid neutral and then continue to clean all the surfaces.

Clean and examine the valve module

⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Make sure not to over-tighten the valve fittings or let liquid spills occur. To prevent valve damage:

1. Keep the instrument clean and dry at all times.
2. Keep all surfaces clean. Use a damp cloth to clean the module surfaces. Dry all surfaces thoroughly. If an acid spill occurs, neutralize the acid with baking soda and then proceed to clean all the surfaces.
3. The rotor seal in the valve module wears with use. Routine replacement of the module is necessary. The main cause of early failure is the presence of abrasive particles in the samples. Filter or centrifuge the sample to prevent scratches on the rotor seal surface.
4. If leaks occur and the above suggestions do not help, contact the manufacturer for a job number and send the unit to the factory immediately to reduce the severity of the damage.


Connect the switching valve

The switching valve is a manual valve that moves the flow in the manifold. Some manifolds that include a column (e.g., cadmium column for nitrate reduction) can have this type of valve.

A flared tubing connection uses flared tubing 0.5 mm ID from the switching valve flare kit. Refer to [Figure 3](#).

Attach the connector to the valve and finger-tighten.

Figure 3 Flared tubing connections

1 Switching valve	3 Metal washer (4x)
2 Valve connector (4x)	4 Flare (4x)


Identify a flow cell leak

If there is a flow cell leak, issues can occur. A leak can be caused by a poorly sealed flare. The flow cell should be examined regularly for leaks and signs of a leak. Do the steps that follow when a leak is found.

1. Remove the flow cell immediately to keep all liquids away from the electronics inside the detector head. Refer to [Figure 4](#).
2. Clean the ports of the flow cell with a cotton swab and deionized water and assemble again.
Note: If the leaks come from the flow cell edges, make sure that there are tightened, good seals at both ports.
3. If air bubbles stay inside the cell, clean the new flow cells. Remove and clean the fiber optic cables.
4. When the manifold is installed, put all the reagent tubing in a 0.25 M sodium hydroxide solution.

- Operate the system with the caustic solution for five minutes.
- Put all the reagent tubing in deionized water to flush out the sodium hydroxide. Operate the system with the deionized water for at least 10 minutes.

Figure 4 Flow cell connections

1 Flow cell	4 Flow cell connector (tubing goes to waste)
2 Flare	5 Flow cell connector (tubing comes from manifold)
3 O-rings (3x)	

Replace the lamp

⚠ DANGER

Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.

⚠ CAUTION


Burn hazard. Turn off power and allow the instrument to cool before this procedure.

NOTICE

Do not touch the lamp bulb when the lamp is installed.

Refer to the specifications in the installation manual for the lamp specifications and refer to [Figure 5](#) for the lamp replacement.

Figure 5 Lamp replacement

Replace the fuse


▲ DANGER	
	Fire hazard. Use the same type and current rating to replace fuses.
▲ DANGER	
	Electrocution hazard. Always remove power to the instrument before making electrical connections.

Refer to the specifications in the installation manual for the user-replaceable fuse specifications.

Replace the system fuse

To replace the fuse, refer to [Figure 6](#).


Figure 6 System fuse replacement

Replace the pump fuse

To replace the fuse, refer to [Figure 7](#).

Figure 7 Pump fuse replacement

Troubleshooting

Problem	Possible cause	Solution
The flow moves backwards.	There is a blockage on the manifold.	Refer to Find the flow blockages on page 18.
Air bubbles come out of the heater block, but no air gets in the heater block.	Dissolved gases come out of the solution	• Prepare reagents a day in advance and keep the reagents at room temperature. • Degas all reagents with helium (unless otherwise noted) • Use a longer backpressure loop. Increase in 50 cm increments.
A white precipitate is formed inside the manifold tubing when TKN is flushed.	Order of the reagent tubing removal	Remove the sodium hydroxide and buffer tubing only after salicylate reagent is no longer present on the manifold.

Problem	Possible cause	Solution
Areas and heights are shown, but there is no concentration for samples.	The system has not been calibrated.	Do a calibration and analyze the samples again.
Patterned air spikes	Time to valve is not correct	If the air spikes occur with every peak in the same location, look at the time for the sample to get to the valve and the load time.
Random air spikes	Dissolved gases come out of the solution	• Prepare reagents a day in advance and keep the reagents at room temperature. • Degas all reagents with helium (unless otherwise noted). • Use a longer backpressure loop. Increase in 50 cm increments. • See if the air is pulled in the reagent flow intermittently. • Make sure that there is sufficient reagent volumes.
	Transmission tubing seal to nipple is not complete	
	Bad reagent	
Sampler initialization error	Sampler fails to initialize	• Under the "Configuration" tab in the software, select "Autosampler". Select "Initialize Autosampler". • Make sure that the power to the sampler is on. • Make sure that sampler cable is not loose. • Make sure that the latest software is installed.
The computer shows that an application error or a runtime error has occurred.	More than one program operates at the same time.	Select Ignore when asked. Close Windows™ and start the application again.
The standards fail	Multiple causes	• Make sure that the check standards are made and added correctly. • Make sure that the peaks are integrated correctly. • Make sure that the reagent preparation and placement are correct. • Make sure that the time to the valve is correct.
Calibration failure	Peak integration is not correct	
	Standard preparation is not correct	
	Placement of sampler is not correct	
	Flow/chemistry problem	
Sampler time out	Problem with sampler arm movement	• Make sure that the sampler arm moves freely. • Under the "Configuration" tab in the software, select "Autosampler". Select "Initialize Autosampler". • Look for loose or unplugged sampler cable.
	Communication error with the system	
The heater element does not get warm.	Heater temperature is not set correctly.	• Make sure that the heater temperature has been set for that channel in Run Properties>Analytes. Make sure that the set point shows correctly in the Omnion channel data display. • Connect the thermocouple • Contact technical support
	Thermocouple is loose or not connected	
	Controller does not work	

Problem	Possible cause	Solution
The baseline is not at zero.	For most chemistries, the baseline will not be at zero.	- Examine the baseline in the Lachat method support data. The baseline should be less than one volt of the support data in the method. If the baseline increases and/or if the blank peaks are negative, the carrier is probably contaminated. - Remove the filter and flow cell while the baseline is shown. If the baseline changes significantly when either component is removed, the issue can be the component or the installation of the component.
The probe continues to hit the cups.	Wrong sample rack Sample racks are not installed correctly. Probe might be bent.	- Make sure that the correct type of rack is selected in the software. - Make sure that the rack is attached correctly to the sampler - Make sure that the probe is not bent. - Make sure that the probe is aligned correctly with the wash bath.
Peaks are smaller than the ones received in the past.	Bad reagent Sample loop is not correct Reagent order is not correct Timing is not correct Interference filter is not correct	- Make sure that the standard preparation, reagent preparation and shelf life are correct. - Make sure that the sample loop and the correct interference filter are located in the correct location. - Make sure that the transmission tubes are in the correct containers. - Make sure that the time for sample to valve is good.
The ammonia peaks are very small.	The pH of the chemistry is not correct The hypochlorite is expired or not made correctly The heater is not set to the correct temperature	- Make sure that the pH at the waste line is 13. - Do the QuikChem method for preparation. Make a fresh solution. - Make sure that the heater is set at 60 °C (140 °F) or other temperature given in the QuikChem method.
Stains in the tubing (after the heater) when ammonia is used	The pH is not correct Caused by the use of household bleach	- Make sure that there is pH at the exit of the flow cell. - Use a reagent-grade 4–6% hypochlorite solution.
The QuikChem method shows support data that cannot be made on the system.	Standard and/or reagent preparation is not correct Manifold configuration is not correct Interference filter and/or sample loop is not correct Flow problems	- Make sure that the formulas in the QuikChem method are done correctly and that the standards have been prepared correctly. - Make sure that the manifold, the interference filter wavelength and the sample loop are correct. - Look for flow problems.
Flow cell leaks at fittings or around glass	Worn out O-rings and flares	- Change the flow cell flares and O-rings. Use the pump to move water through the flow cell to make sure that the leak has stopped. - Clean the leaked fluids off the detector and the fiber optic cables.

Do a dye test

Use a dye test is to find a problem with the reagent flow or hardware setup. If nothing is found, the problem can be chemistry related. Use an interference filter of 500 nm or larger wavelength and a universal dye. Refer to [Make a universal dye](#) on page 17.

There are two important dye tests that can be done. The first test is used to find flow problems. The second test can help find a valve problem or partial clogs that could not be found with the first test. The latter test can also be used to find out whether the sample loop is filled correctly.

Make a universal dye

1. In a one liter container, put 7.5 mL (0.25 fl. oz) of each food dye: green, yellow, red and blue.
2. Fill the container with deionized water and shake it several times.

Note: Universal dye is available from the manufacturer (50959).

Identify flow problems with a dye test

This test will find blockages in the reagent tubing or air leaks in the manifold.

1. Put each transmission tube (includes the carrier and the sample probe) in the universal dye one at a time. Make sure that the dye flows through the tubing. If there is a clog, refer to [Find the flow blockages](#) on page 18.
2. Look for transmission tubes that pulsate. This might be caused by a partial clog. Make sure that the liquid from the waste tubing falls in the waste container and that the waste lines are not below the level of the liquid in the waste container.
Note: The term pulsate refers to the dye that moves forward and backward constantly, but still moves toward the manifold. The tubing can also pulse.
3. If the dye suddenly stops and flows back to the container, look for a clog somewhere in that reagent stream. The path of the solution from the reagent container to the waste container must be located after the flow cell.
4. Look at the transmission tubing and tubing on the manifold to see if there are air bubbles. If yes, look at the connections prior to the location where the bubbles are visible. Look for crimped tubing where the connection is made. There must be two O-rings at every tee or union connector.
5. Look at the flow of the dye to see the correct path of that particular reagent. This test can also show connections that are not correct on the manifold or at the valve. A particular tubing can be in a universal dye solution and the others in deionized water.

Identify valve problems with a dye test

Make a series of dye standards that can be operate with the method in use. Make two injections for each standard and use a first-order-polynomial fit.

Another option is to make up four dye standards and write a new method with exactly the same parameters for the method in use. The concentrations for the standards could be a sequence of four numbers in descending order. Make sure that each standard is half of the preceding one (e.g., 20, 10, 5, 2.5). Refer to [Table 2](#).

1. Put all reagent tubing in deionized water. The manifold for this test agrees with the diagram of the method that is operated. Make sure that the pump operates at 35 and the sample loop is approximately 20 cm of the Teflon 0.8 mm ID.
2. Put the dye standards in the standard vials and start a calibration operation. Monitor the dye as it fills the sample loop when the valve is in the load state. Monitor the dye as it flows to the manifold when the valve is in the inject state.
3. The sample loop should be filled with dye when the valve moves from the load to the inject state. If the sample loop is not filled completely, adjust the Valve Timing.
4. The calibration will pass if the standards were made correctly. If poor precision occurs, do a valve test (refer to [Identify injection valve faults](#) on page 18) and look for clogs (refer to [Find the flow blockages](#) on page 18). A poor precision is specified as a %RSD of more than 10%.
5. If the calibration is correct, the hardware and manifold are in good condition. To make sure that the chemistry is correct, look at the standards and the full method procedure again. Make sure

that the manifold agrees with the diagram in the QuikChem method. Look at the color of the reagent pump tubes, the wavelength of the interference filter, the sample loop length, the heater temperature if it is a heated chemistry, etc.

Table 2 Example of dye standards

Standard	Universal dye (mL)	DI water (mL)
A	100	100
B	50	150
C	25	175
D	12.5	187.5

Troubleshooting based on peak shapes

For screen captures, detailed information and solutions, refer to the expanded user manual available on the manufacturer's website.

Peak shape interpretation	Common causes
Refractive index (Lachat uses this effect in the methods for sea water and brackish water.)	Overacidified samples; wrong carrier
pH effect (such as a dip in the center of the peaks)	Buffer solution with different pH; sample pH, wrong pump tubes; clogs or partial clogs; reagent preparation
Valve artifact (noise that occurs in the baseline when the valve goes from load to inject)	Valve; sample loop; flow restrictions in manifold; waste tubing put in water or does not flow freely
Carryover	Wrong cycle period; wrong sample loop; wrong standard; rinse station issue
Flat-topped peaks	Wrong standards; reagents; wrong sample loop; manifold; overrange samples; wrong interference filter; dirty optics or other detector issue
No peaks	Valve timing settings; reagents; manifold setup; detector; pump power set to off
Noisy baseline	Pump tubes; blockages; injection valve; flow cell; pump; detector; light bulb; reagents; bubble or particle caught in the flow cell; dirty optics or other detector issue
Drifting baseline	Old chemical or reagent; flow problem; contamination; leaking flow cell
All negative peaks or all positive peaks	Carrier contamination; old carrier solution; inverse chemistry

Identify injection valve faults

There are two main problems that can occur with the valve: clogs and hardware failure. To find out if there are blockages in the injection valve, refer to [Find the flow blockages](#) on page 18. Make sure that the cables are connected correctly if the valve does not turn. A loose valve cable can also cause problems. Make sure that the valve is connected correctly. If a motor or electronic failure occurs, replace the valve or send the valve to the manufacturer for repair.

Find the flow blockages

1. If the pump makes an unusual noise, lift each tubing to see if air is aspirated. Examine the carrier tubing to see if there are leaks. If yes, the carrier tubing is the suspect tubing.
2. Examine the nipple by the glass weight. Make sure that there is no blockage to the aperture.

3. Disconnect the carrier tubing fitting that is attached to the carrier pump tubing adapter, which is found between the pump tubing and the valve.
 - a. If the flow problem is still there when the carrier tubing is disconnected, find the clog or crimp location before the disconnected point.
 - b. If the flow moves forward, there are no total clogs up to that point. Put the carrier connection back in place. The flow problem will be back. Continue this process to find the clog or more blockages down stream.
4. Loosen the connector at port 2 of the valve.
 - a. If the flow problem goes away, there are no clogs up to that point before the valve. Connect port 2 and continue to look for the clog.
 - b. If the problem does not go away after port 2 is disconnected, there is a clog between the pump tubing adapter and the valve connector.
5. Loosen the valve connector at port 3.
 - a. If port 3 does not have a leak, there is an internal clog in the valve or a blockage in the sample loop.
 - b. If port 3 has a leak, look at each connection on the manifold with this process. To remove clogs, clean the connectors (tees, unions, etc.) or replace the Teflon tubing.

Sometimes a dye test is necessary to find the a clog or partial clog. Refer to [Do a dye test](#) on page 17.

Identify a detector fault

It is important not to connect a detector to the system unit when the power is set to on. To find out if a detector is at fault:

1. Operate a method. It can be a dry operation.
2. While the method operates, pull the flow cell out to see if there is a response on the diagram. If the detector does not give a result, contact technical support.
Note: *If the detector is unplugged, the baseline will be located at zero.*

Identify the dilutor accuracy

To make sure that the accuracy of the dilutor is correct:

1. Operate five manually-diluted and five autodiluted standards. Compare the results to find out if the problem is the dilutor or elsewhere.
2. Operate a dye test. Use the dilutor to make a series of standards from the stock universal dye. Use one manifold. Make sure that the concentrations entered in the sample tab agree with the dilutions that are used.
3. Make sure that the filter wavelength is 480 to 660 nm and the sample loop is 20 to 25 cm. Make sure that the calibration curve received is linear and the residuals are less than 1 percent.
4. If the dilutor leaks, contact technical support.
Note: *Make sure that the high standard peak amplitude is less than 5 V. If it is higher, change the stock universal dye 1:2 with deionized water. Otherwise, the calibration will not be linear. Make sure that the wavelength of the filter used is less than 660 nm. Peaks are not detected with an 880-nm filter.*

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)