
DOC012.97.90328

SIP 10

04/2026, Edition 4

User Manual
Manual del usuario
Manuel d'utilisation
Manual do Usuário

English 3

Español 15

Français 27

Português 39

Figures ■ Figuras ■ Figures ■ Figuras 51

Table of Contents

1	Specifications	on page 3	5	Maintenance	on page 11
2	General information	on page 3	6	Troubleshooting	on page 13
3	Installation	on page 6	7	Parts and accessories	on page 13
4	Operation	on page 8			

Section 1 Specifications

Specifications are subject to change without notice.

EN

Specification	Details
Dimensions (W x H x D)	120 × 85 × 200 mm (4.72 x 3.35 x 7.87 inches)
Enclosure	IP30
Weight	0.5 kg (11 lb)
Pollution degree	2
Overvoltage category	I
Protection class	Class III
Cell	Plastic: Dual Path Length Pour-Thru Cell, 1 inch/1 cm Quartz glass: Pour-Thru Cell, 1 cm
Wavelength range	Plastic: 340 to 900 nm Quartz glass: 190 to 900 nm
Rinsing volume	Minimum of 25 mL
Flow speed	1 mL/second
Interface	USB type A
Power requirements	USB cable, 530 mA, 5 V
Operating temperature	10 to 40 °C (50 to 104 °F), 80% relative humidity, non-condensing
Storage temperature	−10 to 60 °C (14 to 140 °F), 85% relative humidity, non-condensing
Certifications	CE, UKCA, CMIM
Warranty	1 year (EU: 2 years)

Section 2 General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual unless otherwise required by applicable law or contract between the parties. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

2.1 Safety information

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

If the equipment is used in a manner that is not specified by the manufacturer, the protection provided by the equipment may be impaired. Do not use or install this equipment in any manner other than that specified in this manual.

2.2 Use of hazard information

EN	⚠ DANGER Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.
	⚠ WARNING Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.
	⚠ CAUTION Indicates a potentially hazardous situation that may result in minor or moderate injury.
	NOTICE Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

2.3 Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

2.4 Certification

Canadian Radio Interference-Causing Equipment Regulation, ICES-003, Class B:

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FCC Part 15, Class "B" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Move the equipment away from the device receiving the interference.
2. Reposition the receiving antenna for the device receiving the interference.
3. Try combinations of the above.

2.5 Chemical and biological safety

Chemical or biological hazards. If this instrument is used to monitor a treatment process and/or chemical feed system for which there are regulatory limits and monitoring requirements related to public health, public safety, food or beverage manufacture or processing, it is the responsibility of the user of this instrument to know and abide by any applicable regulation and to have sufficient and appropriate mechanisms in place for compliance with applicable regulations in the event of malfunction of the instrument.

⚠ DANGER
Fire hazard. This product is not designed for use with flammable liquids.

- Do not use the instrument in samples that could cause a biological hazard.
- Observe all cautionary information printed on the original solution containers and safety data sheets prior to their use.
- Dispose of all consumed solutions in accordance with the local and national regulations and laws.
- Select the type of protective equipment suitable to the concentration and quantity of the dangerous material being used.

2.6 Product overview

The SIP 10 Sipper module is an optional accessory for the DR3900, DR4900 and DR6000 spectrophotometer. The measurement accuracy is improved by working with the SIP 10. The same optical characteristics apply to both the zero compensation and the measurements. All test solutions flow through a cell and errors due to optical differences between different cells are eliminated. A peristaltic pump pumps a constant amount of fluid through the Pour-Thru Cell.

2.7 Product components

Make sure that all components have been received. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

DR3900 and DR4900

- SIP 10 Sipper module, pump tubing pre-assembled
- Dual Path Length Pour-Thru Cell, 1 inch/1 cm
- USB connection cable

- Tygon® tubing, 1.70 m (5.6 ft)
- Stainless steel tube adapter

DR6000

- SIP 10 Sipper module, pump tubing pre-assembled
- Pour-Thru Cell according to order variant:
 - Dual Path Length Pour-Thru Cell, 1 inch/1 cm, plastic
 - Pour-Thru Cell, 1 cm, quartz glass
- USB connection cable
- Tube set consisting of
 - 2× Tygon tubes with fitting
 - Tygon tube
 - 2× Pharmed® tubes with fitting
 - Rubber tube guide
- Stainless steel tube adapter
- Cell compartment cover
- Tray for storage

Section 3 Installation

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

⚠ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

NOTICE

Make sure that tubes are not bent. Check the fill level of the waste vessel regularly.

NOTICE

Connect the device directly to the photometer via the USB cable. Do not use a USB hub connection.

3.1 DR3900 or DR4900 installation

1. Position the SIP 10 to the left of the photometer. Refer to [Figure 1](#) on page 51 and [Figure 2](#) on page 52.
2. Clean the Pour-Thru Cell with a lint-free cloth.
3. Place the Pour-Thru Cell into the cell compartment of the photometer according to the required path length.
4. Cut a piece of tubing no longer than 40 cm (1.3 ft) for the feed tube.
5. Connect one end of the feed tube to the IN connection of the Pour- Thru Cell.
6. Connect the other end of the feed tube to the stainless steel tube adapter and place this in the zero-solution/sample vessel.
7. Cut a piece of tubing approx. 30 cm (1 ft) long for the connection tube.
8. Connect the connection tube to the OUT connection of the Pour- Thru Cell and the input connection of the peristaltic pump.

9. Use the remaining tubing as a discharge tube, and connect one end to the output connection of the peristaltic pump.
10. Place the other end of the discharge tube in a suitable waste vessel.
11. Connect the SIP 10 USB cable to a USB port on the photometer.

3.2 DR6000 installation

1. Position the SIP 10 on the storage tray to the left of the photometer. Refer to [Figure 3](#) on page 53.
2. Open the cell compartment cover up to the stop.
3. Push the rubber tube guide into the groove in the top left of the cell compartment.
4. Push the two Pharmed tubes with the free ends onto both connections of the Pour-Thru Cell.
5. Clean the Pour-Thru Cell with a lint-free cloth.
6. Place the Pour-Thru Cell into the cell compartment of the photometer according to the required path length.
7. Push the Pharmed tube with the connection to the IN-connection of the Pour-Thru Cell into the bottom guide in the tube guide.
8. Push the Pharmed tube with the connection to the OUT-connection of the Pour-Thru Cell into the top guide in the tube guide.
Note: *Make sure that the tubes are sitting securely in the guide so that they do not become squashed.*
 The fitting ends of two Pharmed tubes are now protruding from the side of the photometer.
9. Install the cell compartment cover, starting at the bottom right corner.
10. Close the cover.
11. Connect the fitting end of the bottom Pharmed tube to the fitting side of a Tygon tube.
 This is the feed tube.
12. Connect the other end of the feed tube to the stainless steel tube adapter and place this in the zero-solution/sample vessel.
13. Connect the fitting end of the top Pharmed tube to the fitting side of the second Tygon tube.
 This is the connection tube.
14. Connect the other end of this connection tube to the input connection of the peristaltic pump.
15. Connect one end of the Tygon tube (without fitting) to the output connection of the peristaltic pump.
 This is the discharge tube.
16. Place the other end of the discharge tube in a suitable waste vessel.
17. Connect the SIP 10 USB cable to a USB port on the photometer.
Note: *Shorten the feed tube and connection tube as much as possible in order to optimize the pump volume of the system.*

Section 4 Operation

4.1 SIP 10 setup for DR3900 or DR6000

When a SIP 10 module is connected to the photometer, an additional **Sipper** button is shown in the toolbar.

1. Press **Sipper** on the toolbar. A Sipper options menu opens.
2. Press **On** to modify the settings for a Sipper cycle.

EN

Option	Description
Sip time	The sip time determines the sample volume that is pumped through the cell. In order to make sure that all of the liquid is replaced in the Dual Path Length Pour-Thru Cell, 1 inch/1 cm, the minimum time is 25 seconds. Default: 25 seconds. The flow rate is 1 mL/sec. Note: If you use other Pour-Thru Cells, the SIP time must be determined individually. One method to determine this is to use two solutions with different concentrations and measure the absorbance continuously.
Settle time	The settle time defines the duration of the standstill phase between the pump process and measurement process. During this time, air bubbles that have formed during the pump process can escape and sample turbulence can cease. Default setting: 5 seconds
Purge time	The purge time determines the volume of rinse agent that is pumped through the cell after each measurement. This additional purging can be deactivated. Default setting: off
Purge start	The start of the purge time can be triggered automatically after the measurement process, or manually. To start the cycle manually, press Purge to trigger the purge cycle. Default setting: auto Note: The Purge start button is only active when the purge time setting is On .

3. Configure the required settings and press **OK** to confirm.

4.1.1 Execution of a stored program

For the tests listed in [Table 1](#) on page 8, an additional purge cycle of the Pour-Thru Cell with deionized water is required between samples.

Table 1 Tests for which purging with deionized water is required

Aluminum, Aluminon	Chlorine dioxide, LR	Cobalt, PAN
Copper, Porphyrin	Hardness, Calmagite	Manganese, LR, PAN
Nickel, PAN	Nitrate, MR	Nitrate, HR

The Dual Path Length Pour-Thru Cell, 1 inch/1 cm can also be used for the Nessler Method for nitrogen and TKN if cleaned correctly. For this, pour a few sodium thiosulfate crystals into the cell to clean it. Rinse out the crystals with deionized water. The use of the Pour-Thru Cell is not possible for the tests listed in [Table 1](#) on page 8, since chemical problems or other complications can occur. Use the sample cell specified in the procedural instructions for these tests.

Table 2 Tests for which the Pour-Thru Cell is not suitable

Aluminum ECR	Arsenic	Barium	Boron, Carmine
Cyanuric acid	Fluoride	Formaldehyde	Lead, LeadTrak
Mercury	Nickel, Heptoxime	Nitrite, HR	PCB
Phenols	Potassium	Selenium	Silver
Suspended solids	Sulfate	TPH	Volatile acids
Zinc		Surfactants, Anionic (Detergents)	

The Pour-Thru Cell can also be used for other tests. However, please note that the Pour-Thru Cell requires a minimum sample volume of 25 mL to make sure that the previous solution is fully rinsed out of the cell. Increase the reagent and sample volumes accordingly.

1. Access the required test via **Stored programs**.

2. Place the feed tube in the zero solution and press **Zero**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.

Note: If the **Purge start: Auto** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

3. Place the feed tube in the sample solution and press **Read**.

The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.

Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.

Note: If the **Purge start: Auto** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

Note: Either deionized water or the next sample can be used as a rinse agent.

The measurement result is shown on the display.

4.1.2 Execution of a single or multi-wavelength measurement

1. Press **Single wavelength** or **Multi wavelength**.

2. Press **Options** to adjust the settings within single/multiwavelength mode.

3. Place the feed tube in the zero solution and press **Zero**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.

Note: If the **Purge start: Auto** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

4. Place the feed tube in the sample solution and press **Read**.

The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.

Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.

Note: If the **Purge start: Auto** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

Note: Either deionized water or the next sample can be used as a rinse agent.

The measurement result is shown on the display.

4.1.3 Execution of a time course measurement

1. Press **Time course**.

2. Press **Options** to adjust the settings within time course mode.

3. Place the feed tube in the zero solution and press **Zero**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: If the **Purge start: Auto** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

4. Place the feed tube in the sample solution and press **Read**.

The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.

The measurement is taken in accordance with the time course settings. The **Zero** and **Read** buttons change to **Mark** and **Stop** during the time course. To abort the measurement process, press **Stop**.

Once the entire time has passed, the measurement is shown as a curve on the display.

Note: If the **Purge start: Auto** setting has been activated in the **Sipper options** menu, the rinse agent is pumped in automatically after measurement.

4.2 SIP 10 setup for DR4900

1. Push , then select **Settings > Measurement settings**.
2. Scroll down to the **Sipper** options.
3. Set the menu to on when the sipper module is connected.
4. Configure the applicable settings:

Option	Description
Sipper time (in seconds)	Sets the sip time to select the sample volume that is pushed through the sample cell. Default: 8 seconds
Settle time (in seconds)	Sets the time between the pump procedure and measurement procedure (standstill interval). Default: 5 seconds
Purge	Set to on to select rinse options.
Purge time (in seconds)	Sets the time period during which the detergent is pushed through the sample cell after each measurement. Default: Off
Purge start	Select Automatic (default) to start the purge procedure after each measurement. Select Manual to manually start a purge procedure when necessary. Note: When Manual is selected, a Purge button is shown.
Purge test	Push Start to start the test and rinse procedure.

4.2.1 Execution of a Hach program

1. Push , then select **Sample analysis > Hach programs**.
2. Select the applicable procedure from the list and push **Start**.
3. Place the feed tube in the zero solution and press **Zero**.
The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.
Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.
Note: If the **Purge start > Automatic** setting has been activated in the **Sipper options** menu, the rinse agent is pumped in automatically after measurement.
The zero measurement is shown on the display.
4. Place the feed tube in the sample solution and press **Read**.
The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.
Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.
Note: If the **Purge start > Automatic** setting has been activated in the **Sipper options** menu, the rinse agent is pumped in automatically after measurement.
Note: Either deionized water or the next sample can be used as a rinse agent.
The measured values show on the display. The measured values are automatically saved to the data log if **Auto save** is enabled.

4.2.2 Execution of a single or multi-wavelength measurement

1. Push , then select **Sample analysis > Single wavelength** or **Multi wavelength**.
2. Place the feed tube in the zero solution and press **Zero**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.

Note: If the **Purge start > Automatic** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

3. Place the feed tube in the sample solution and press **Read**.

The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.

Note: The remaining time is shown on the display. Press **Cancel** to stop the operation.

Note: If the **Purge start > Automatic** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

Note: Either deionized water or the next sample can be used as a rinse agent.

The measurement result is shown on the display.

4.2.3 Execution of a time course measurement

1. Push , then select **Sample analysis > Time course**.

2. Place the feed tube in the zero solution and press **Zero**.

The zero solution is pumped in according to the selected SIP time. Once the settle time has passed, the zero measurement is taken.

Note: If the **Purge start > Automatic** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

The zero measurement is shown on the display.

3. Place the feed tube in the sample solution and press **Read**.

The sample solution is pumped in over the selected SIP time. Once the settle time has passed, the measurement is taken.

Once the entire time has passed, the measurement is shown as a curve on the display.

Note: If the **Purge start > Automatic** setting has been activated in the Sipper options menu, the rinse agent is pumped in automatically after measurement.

Section 5 Maintenance

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

⚠ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

⚠ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

5.1 Clean the instrument

⚠ DANGER

Electrocution hazard. Remove power from the instrument before this procedure is started.

NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Clean the exterior of the instrument with a moist cloth and a mild soap solution.

5.2 Tray for storage (for DR 6000 only)

If you no longer require the SIP 10 for current measurements, you can store all of the accessories in the tray:

- Cell compartment cover and USB connection cable to the right of the SIP 10
- Fitted tubes and rubber tube guide behind the SIP 10
- Pour-Thru Cell in the SIP 10

5.3 Pour-Thru Cell

Purge the Pour-Thru Cell before and after every series of tests with deionized water.

If the Pour-Thru Cell is heavily contaminated, repeat the purge cycle several times.

Before the Pour-Thru Cell is put into storage, perform the purge cycle several times with air to remove any residual fluid. Hold the Pour-Thru Cell with the OUT connection slanted downwards during this procedure.

Check the Pour-Thru Cell viewing windows. If they are dirty or fogged up, place the Pour-Thru Cell in a soap solution or in dilute acid. Then rinse the Pour-Thru Cell thoroughly with deionized water.

5.4 Tube cleaning

⚠ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Clean the tubes with deionized water after each series of measurements.

The tubes are exposed to chemicals and must be replaced at least every 12 months.

5.5 Replace the pump tubing

The pump tubing is exposed to mechanical loads and chemicals and must be replaced at least every 12 months.

1. Perform the purge cycle several times with air to remove any residual fluid from the system.
2. Disconnect the USB cable from the photometer.
3. Remove the connection tube and discharge tube from the input and output connections on the peristaltic pump.
4. Tilt back the lever behind the peristaltic pump. Refer to [Figure 4](#) on page 55, step 1. The peristaltic pump cover will move back.
5. Lift up the peristaltic pump cover and dispose of the cover, pump tubing and connections. Refer to [Figure 4](#) on page 55, step 2.
6. Position the new cover with pre-assembled pump tubing and connections on the peristaltic pump.

7. Push down the cover on the peristaltic pump and tilt up the lever behind the peristaltic pump. Refer to [Figure 4](#) on page 55, step 3 and 4.
The peristaltic pump cover will move forward.
8. Connect the connection tube and the discharge tube to the input and output connections on the peristaltic pump. Refer to [Figure 1](#) on page 51, [Figure 2](#) on page 52 and [Figure 3](#) on page 53.
9. Connect the SIP 10 USB cable to a USB port on the photometer.

Section 6 Troubleshooting

Potential SIP 10 errors are displayed by the photometer. Refer to [Table 3](#) on page 13.

Table 3 Error messages

Error	Possible cause	Solution
Sipper module not connected. Please check cable.	Connection between SIP 10 and photometer is disrupted.	Check the USB cable. The cable length must not exceed 1 m (3.3 ft) and the cable must be directly connected to the photometer. Remove other devices that are connected where necessary.
Please check sipper and tube.	Pump tubing is not correctly inserted.	Loosen the cover and position again. Where required, run the pump briefly before the lever is moved to make sure that the tube is correctly positioned around the rollers.

Section 7 Parts and accessories

⚠ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
SIP 10 Sipper module set for DR 3900 or DR4900 complete with tube set and 1 inch/1 cm Pour-Thru Cell, EU	LQV157.99.10001
SIP 10 Sipper module set for DR 6000 complete with tray, tube set and 1 inch/1 cm Pour-Thru Cell, EU	LQV157.99.20001
SIP 10 Sipper module set DR 6000 complete with tray, tube set and 1 cm quartz glass Pour-Thru Cell, EU	LQV157.99.30001
SIP 10 Sipper module; including pump tubing, EU	LQV157.99.00001
USB cable, type AB 1 m (3.3 ft)	LZQ104
Dual Path Length Pour-Thru Cell, plastic, 1 inch/1 cm	LZV876
UV Pour-Thru Cell, quartz glass, 1 cm	LZV510
Pump tubing (Lagoprene®) pre-assembled with cover and connections	LZV877
Complete tube set for DR 3900, including Tygon tubing, 1.70 m (5.6 ft), pump tubing (Lagoprene pre-assembled with cover and connections), stainless steel tube adapter	LZV875
Complete tube set for DR 6000 and Dual Path Length Pour-Thru Cell, including Tygon tubing and Pharmed tubing, pump tube (Lagoprene pre-installed with cover and connections), stainless steel tube adapter	LZQ102

Replacement parts (continued)

Description	Item no.
Complete tube set for DR 6000 "drinking water application," including Pharmed tubing, pump tube (Lagoprene preinstalled with cover and connections), stainless steel tube adapter	LZQ100
Tray for storage	HTT057
Cell compartment cover	LZQ105
Information on Tygon tubes, EC no. 1907/2006 (REACH)	DOC063.98.90334

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


EN	1 SIP 10 Sipper module on the tray for storage	8 Rubber tube guide
	2 Discharge tube: Tygon tube (without fitting)	9 Pour-Thru Cell
	3 Peristaltic pump	10 Feed tube: Pharmed tube with fitting
	4 USB cable	11 Feed tube: Tygon tube with fitting
	5 Connection tube: Tygon tube with fitting	12 Photometer
	6 Connection tube: Pharmed tube with fitting	13 Stainless steel tube adapter
	7 Cell compartment cover	
ES	1 Módulo Sipper SIP 10 de la bandeja para almacenamiento	8 Guía de manguera de caucho
	2 Manguera de descarga: manguera Tygon (sin conexión)	9 Pour-Thru Cell
	3 Bomba peristáltica	10 Manguera de alimentación: manguera de Pharmed sin conexión
	4 Cable USB	11 Manguera de alimentación: manguera Tygon con conexión
	5 Manguera de conexión: manguera Tygon con conexión	12 Fotómetro
	6 Tubo de conexión: Tubo Pharmed con racor	13 Adaptador para manguera de acero inoxidable
	7 Cubierta del compartimento de cubetas	

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)