Sonda 15 50 T

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7

50 15 T probe

Additional information is available on the manufacturer's website.

Caution
⚠ *Personal injury hazard. Glass components can break. Handle with care to prevent cuts.*
⚠ *Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.*

Introduction
This pH refillable probe with temperature sensor is used in alkaline water samples, samples with proteins and samples with sudden temperature change.

Technical information

Range	0 to 14 pH
Junction	Ceramic (3)
Electrolyte	KCl, 3 M+glycerine
Reference	Ag ion barrier cartridge
Temperature	Pt 1000
Minimum depth	18 mm (0.7 in.)
Operating temp	-10 to 100 °C (14 to 212 °F)
Limits	Viscous solutions, distilled water, samples with colloids

1 Connect the probe

2 Calibration
⚠ *Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.*

3 Measurement

Maintenance and storage
Cleaning: soak the probe in a cleaning solution. Rinse with deionized water.

Contamination	Cleaning agent
Proteins	Pepsin cleaning solution
Grease, oils, fats	Electrode cleaning solution
Limescale	0.1 N HCl solution

Hydration: soak the dry probe in a dilute HCl solution for several hours.
Storage: put the probe in the protector with KCl 3 M and glycerol solution.

Troubleshooting
Be sure that the samples are suitable for the application.

Warranty
6 months for manufacturing faults only.
Damage from use is not covered.

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