**LANGE** 

DOC022.98.90466

# TitraLab pH & Alkalinity

08/2013, Edition 1A

|                                |    |
|--------------------------------|----|
| Basic User Manual              | en |
| Basis Bedienungsanleitung      | de |
| Manuale utente di base         | it |
| Manuel d'utilisation de base   | fr |
| Manual básico del usuario      | es |
| Manual do utilizador básico    | pt |
| Basisgebruikershandleiding     | nl |
| Grundlæggende brugervejledning | da |
| Podstawowa instrukcja obsługi  | pl |
| Grundläggande bruksanvisning   | sv |

|                  |     |
|------------------|-----|
| English .....    | 0   |
| Deutsch .....    | 14  |
| Italiano .....   | 30  |
| Français .....   | 46  |
| Español .....    | 62  |
| Português .....  | 78  |
| Nederlands ..... | 94  |
| Dansk .....      | 110 |
| Polski .....     | 124 |
| Svenska .....    | 140 |

# Specifications

Specifications are subject to change without notice.

| Performance specifications |                     |                                                                                                                                                                               |
|----------------------------|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Measuring principle        |                     | Acid–base titration, pH: potentiometric determination                                                                                                                         |
| Application area           |                     | Determination of pH and Alkalinity in water                                                                                                                                   |
| Measuring range            |                     | Min. 0.1 mL or 10 mg/L CaCO <sub>3</sub> , Recommended max. 10 mL or 1000 mg/L CaCO <sub>3</sub><br>Maximum titration system capacity of 20 mL or 2000 mg/L CaCO <sub>3</sub> |
| Resolution                 |                     | Dependent on the selected units                                                                                                                                               |
| Electrode                  |                     | pH electrode, combined reference electrode, with built-in Pt 1000 sensor                                                                                                      |
| Calibration                |                     | One, two or three-point calibration                                                                                                                                           |
| Buffer recognition         |                     | Automatic: pH 2.00; 4.01; 7.00; 9.21 and 10.90 (at 25 °C)                                                                                                                     |
| Programs                   |                     | pH measurement, 1 titration program, 1 calibration program                                                                                                                    |
| Languages                  |                     | English, Spanish, French and Italian                                                                                                                                          |
| Data logger                |                     | 55 data records                                                                                                                                                               |
| Date/time                  |                     | DD-MM-YYYY                                                                                                                                                                    |
| Display                    |                     | 128 × 64 pt, backlit liquid crystal display                                                                                                                                   |
| Keypad                     |                     | 7 membrane keys, splash-proof, PET, lifetime of >6 million key-presses                                                                                                        |
| Burette                    | Resolution          | 0.001 mL                                                                                                                                                                      |
|                            | Material of syringe | Polypropylene                                                                                                                                                                 |
|                            | Volume              | 10 mL                                                                                                                                                                         |
|                            | Tubes               | PTFE tubes                                                                                                                                                                    |

|                                                                        |                               |                                                                                                                                                        |
|------------------------------------------------------------------------|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Inputs and outputs</b>                                              | <b>Measuring electrode</b>    | BNC connector; maximum voltage $\pm 2.2$ V, maximum current 1 mA. <sup>1</sup>                                                                         |
|                                                                        | <b>Magnetic stirrer</b>       | RCA connector (cinch), maximum voltage 12 V, maximum current 150 mA <sup>2</sup>                                                                       |
|                                                                        | <b>Computer or printer</b>    | RJ11 telephone connector                                                                                                                               |
|                                                                        | <b>Power supply, external</b> | Instrument input: 24 V $\approx$ 45-65 Hz, 1 A<br>External power supply delivered:<br>Input: 100-240 VAC, $\pm 10$ %, 50-60 Hz, Output: 24 VDC, 1.33 A |
|                                                                        | <b>Mouse</b>                  | Mini-DIN connector                                                                                                                                     |
| <b>Dimensions of measuring instrument and environmental conditions</b> |                               |                                                                                                                                                        |
| <b>Enclosure material</b>                                              |                               | ABS and enameled steel                                                                                                                                 |
| <b>Weight</b>                                                          |                               | approximately 4 kg                                                                                                                                     |
| <b>Dimensions</b>                                                      |                               | 130 × 160 × 300 mm                                                                                                                                     |
| <b>Working temperature</b>                                             |                               | +15 to +40 °C (+60 to +104 °F)                                                                                                                         |
| <b>Storage temperature</b>                                             |                               | −10 to +50 °C (+14 to +122 °F)                                                                                                                         |
| <b>Air humidity</b>                                                    |                               | <80 % relative humidity, non-condensing                                                                                                                |
| <b>Altitude requirements</b>                                           |                               | Standard 2000 m (6562 ft) ASL (Above Sea Level)                                                                                                        |
| <b>Certification</b>                                                   |                               |                                                                                                                                                        |
| <b>Certification</b>                                                   |                               | CE                                                                                                                                                     |
| <b>Electrical security</b>                                             |                               | According CE, EN 61010                                                                                                                                 |
| <b>Electromagnetic compatibility</b>                                   |                               | According CE, EN 50081-2 and EN 50082-2                                                                                                                |
| <b>Warranty</b>                                                        |                               |                                                                                                                                                        |
| <b>Warranty</b>                                                        |                               | 2 years                                                                                                                                                |

<sup>1</sup> Only for supplied electrode and BCN cable.

<sup>2</sup> Only for supplied stirrers.

## General information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

## Safety information

Please read the entire user manual carefully before the equipment is unpacked, set up or put into operation. Pay attention to all danger and warning statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

To make sure that the protection provided by this instrument is not impaired, do not use or install this instrument in any manner other than that specified in these operating instructions.

### ! DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

### ! WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

### ! CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

**Note:** Information that supplements points in the main text.

## Warning labels

Read all labels and tags attached to the instrument. Failure to do so may result in personal injury or damage to the instrument.

For symbols applied to the instrument, corresponding warnings are provided in the user manual.

This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.

This symbol may be found on an enclosure or barrier within the product and indicates a risk of electrical shock and/or death by electrocution.

Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems after 12 August 2005. In conformity with European local and national regulations, European electrical equipment users must now return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

**Note:** For return for recycling, please contact the equipment manufacturer or supplier for instructions on how to return end-of-life equipment, manufacturer-supplied electrical accessories, and all auxiliary items for proper disposal.

## Safety information for the correct use of the instrument

- This instrument is intended only for indoor use.
- Do not vibrate or jolt the instrument.
- Do not open the instrument. Only authorized and trained personal is allowed to open the instrument. Follow the following safety information, in addition to any local guidelines in force.
- Guarantee is voided if the instrument is not used in accordance with the guidelines present in this document.

## Chemical and biological safety

### WARNING

Potential danger in the event of contact with chemical/biological materials.

Handling chemical samples, standards and reagents can be dangerous.

Familiarize yourself with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

### WARNING

This instrument has been developed to meet the requirements of volumetric titration applications. It is therefore aimed at experienced users who have the knowledge required to operate the instrument and implement the security instructions enclosed. Remember that this instrument must not, under any circumstances, be used to perform test on living beings.

Normal operation of this instrument may involve the use of hazardous chemicals or biologically harmful samples.

- Observe all cautionary information printed on the original solution containers and safety data sheet prior to their use.
- Dispose of all consumed solutions in accordance with national regulations and laws.
- Select the type of protective equipment suitable to the concentration and quantity of the dangerous material at the respective work place.

## Installation

### DANGER

Only qualified personnel should conduct the tasks described in this section of the manual.

Safety of any system incorporating this instrument is responsibility of the assembler of the system.

### WARNING

Risk of injury to eyes, skin and respiratory

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Make sure that the system is only unpacked, assembled, connected and operated by qualified and/or trained personnel.


## Unpack the instrument

Unpack all supplied parts carefully, as they are highly sensitive in part to shock and impact. Read the user manual prior to installation and proceed exactly as described.

## Scope of delivery

Please check that the order is complete. If anything is missing or damaged, please contact the manufacturer or retailer.

**Figure 1 Scope of delivery**

|          |                                                   |           |                                        |
|----------|---------------------------------------------------|-----------|----------------------------------------|
| <b>1</b> | TitraLab pH & Alkalinity content system           | <b>8</b>  | pH electrode with cable                |
| <b>2</b> | Tube with bottle closure, DIN 45 thread           | <b>9</b>  | O-ring to secure tube                  |
| <b>3</b> | Tube to attach to electrode                       | <b>10</b> | Adapter for power supply               |
| <b>4</b> | Titration stand holder including magnetic stirrer | <b>11</b> | Power supply, external                 |
| <b>5</b> | 10 × titration beakers, PP, 100 mL                | <b>12</b> | Hardware sleeve                        |
| <b>6</b> | 10 × Teflon magnetic stirring bars (20 mm × 6 mm) | <b>13</b> | Syringe with plunger                   |
| <b>7</b> | 10 × white O-rings for burette and pump threads   | <b>14</b> | Empty desiccant cartridge with stopper |

## Operating environment

### WARNING

The instrument may not be used in dangerous environments.  
The manufacturer and its suppliers reject any express or indirect guarantee for use with high-risk activities.

### NOTICE

Protect the instrument from temperature extremes, including heaters, direct sunlight and other heat sources.

Note the following points so that the instrument will operate perfectly and have a long service life.

- Place the instrument firmly on an even surface. Do not push any objects under the instrument.
- The ambient temperature must be +15 to +40 °C (+60 to +104 °F).
- The relative humidity should be less than 80 %; moisture should not condense on the instrument.
- A minimum clearance of 15 cm (5.9 in.) must be ensured above and on all sides of the device; this allows the air to circulate and prevents overheating of the electric parts. The dosage arm of the stepper motor must be able to move freely.
- Do not operate or store the instrument in extremely dusty, humid or wet locations.
- Make sure that no liquids enter the instrument and immediately wipe away any liquids that contact the instrument.
- Protect the instrument against vibrations, direct sunlight, corrosive gases as well as strong magnetic and/or electric fields.
- Only use genuine replacement parts and accessories.
- Install the instrument in a levelled place and on a compact surface to minimize the sound level of the instrument. This will avoid undesired resonances that can increase the sound level.

- Keep the surface of the instrument and all accessories clean and dry at all times. Splashes or spills on the instrument should be cleaned up immediately.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

## Installation

### Installation of the burette

The burette dispenses the titrant: Hydrochloric acid (HCl) 0.1 N (0.1 M) or Sulfuric acid (H<sub>2</sub>SO<sub>4</sub>) 0.1 N (0.05 M).

The burette comprises the syringe, stepper motor and valve. The stepper motor and valve are already installed in the TitraLab; the syringe has to be installed by the user.

### WARNING

Risk of injury to pinch points and cuts

Make sure that the titrator is never used/operated without the syringe cover.

Do not introduce your fingers on the bottom of the cover while the instrument is connected.

Make sure that the system is only assembled by qualified and trained personnel.

### NOTICE

Warning against system damage

Make sure that the syringe is never operated without the attachment sleeve.

Only use genuine replacement parts and accessories.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

**Note:** Make sure that you do not misplace the attachment sleeve for the syringe. It is not provided with a new syringe.

1. Push the syringe piston into the syringe body.  
Slide the attachment sleeve fully onto the syringe so that the narrow end points toward the piston.
2. Remove the syringe cover by unscrewing it.

3. Point the syringe tip upward, place the syringe into the bracket such that the attachment sleeve is above the bracket.
4. Hold the attachment sleeve in the bracket and screw the syringe onto the valve from below. With the narrow side facing down, carefully press the attachment sleeve into the bracket.
5. Withdraw the piston from the syringe and place it onto the screw in the dosage arm. Hold the syringe piston and turn the attachment screw on the dosage arm into the thread in the piston.
6. Fix the syringe cover with six screws.


## ⚠ WARNING

Risk of injury to eyes, skin and respiratory

During the first days of system use, check the screw joint between the syringe and valve to avoid leaks and the potential escape of reagents.

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

**Figure 2 Installation of the burette**

## Installation of the titration stand

The titration stand consists of a holder onto which the electrode, titration cell and magnetic stirrer are attached. The titration is performed in the titration cell.

1. Place the titration stand on the right side (syringe side) of the titrator (observe the routing of the cable).
2. Place a magnetic stirring bar into a titration cell. Place a titration cell over the magnetic stirrer.
3. Place the rubber ring on the stand rod to set the minimum height of the electrode holder.

- Press the button to lock the electrode holder and place the electrode holder on the stand rod. The electrode holder is locked upon release of the button.
- Remove the white cap from the tip of the electrode. Insert the electrode into the electrode holder.

## NOTICE


### Warning against system damage

Make sure that the height setting of the electrode holder protects the electrode against damage from the stirring bar and against contact with the bottom of the titration cell.

Keep the lock button pressed whenever the electrode holder is being moved.

**Note:** Whenever the electrode is taken from the titration stand or the titration cell is removed from the titration stand, always place a container beneath the electrode. Otherwise, drips from the sample or chemicals could contaminate the work area.

**Figure 3 Installation of the titration stand**

## Electrical connection of the magnetic stirrer

- Connect the cable of the magnetic stirrer to the **STIRRER** port on the back of the TitraLab. The power is supplied via the TitraLab.

## Installation of tubes

## ⚠ WARNING

Risk of injury to eyes, skin and respiratory

During the first days of system use, check the screw joint between the syringe and valve to avoid leaks and the potential escape of reagents.

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Make sure that the tubes are never kinked.

Regularly check the system as well as all tubes and connections for leaks, condition and cleanliness.


To prevent twisting of the tubes, you should always connect flange tube connections before those with grommets.

The system is delivered with two tubes; these are supplied pre-installed on the TitraLab and only need to be connected to the relevant titrant bottle and to the electrode.

## Removal of protective stickers (Uninstalled tubes configuration)

- Remove the protective stickers on the burette valve.
- Make sure that the pre-installed white O-rings are not removed.
- Check that the O-rings are in a horizontal position inside the thread.

Figure 4 Position of protective stickers


### Connection of feed tube

The system is delivered with a tube with bottle closure (DIN 45 thread). This feed tube is supplied pre-installed on the TitraLab and only needs to be connected to the titrant bottle.

1. Connect the tube from the IN socket of the burette valve to the titrant bottle (HCl or H<sub>2</sub>SO<sub>4</sub>). Screw on the bottle closure to do this.

Figure 5 Connection of feed tube

### Connection of the electrode


1. Connect the electrode cable to the **ELECTRODE** port on the back of the TitraLab.  
The power is supplied via the TitraLab.

### Connection of the tube to the electrode

The sensor holder has a tube guide.

1. Connect the tube on the OUT socket of the burette valve.
2. Place the electrode in the sensor holder.
3. Insert the tube into the cannula of the sensor holder. Push the free end of the tube into the tube guide of the sensor holder until the end of the tube is flush with the tip of the electrode (1).
4. Attach the tube to the electrode with a rubber ring (2).

**Figure 6 Connection of the electrode**

## Power supply, external

### **⚠ DANGER**

Risk of electric shock

Make sure the AC supply voltage is suitable for the power supply and the mains plug style is suitable for the type of outlet socket.

Only use the original power supply and plug connectors.


### **⚠ DANGER**

Fire and electrical shock hazard.

Take care, when using chemical products and liquids, to avoid dropping them on the power supply or mains connections.

The power supply is provided with different plug connectors. The correct plug connector must be attached to the power supply before first use.

**Figure 7 IPlug connector attachment**

**Note:** To change the plug connector, press in the catch and slide the plug connector upward.

## Keypad

### Navigation key UP

- Used to select between different options on the display.
- Increases numeric values.
- Shows the previous data record in the data log menu.
- Empties the syringe through the discharge tube during manual activation.

### Navigation key DOWN

- Used to select between different options on the display.
- Reduces numeric values.
- Shows the next data record in the data log menu.

- Fills the syringe through the feed tube during manual activation.

## ON/OFF key

### Navigation key LEFT

- Moves the cursor backward for entry of numeric values.
- Scrolls up in menus.
- Empties the syringe through the feed tube during manual activation.

### Navigation key RIGHT

- Moves the cursor forward for entry of numeric values.
- Scrolls down in menus.

## RETURN key

- Scrolls back in menus.
- Cancels a manual control.
- Cancels a titration.

## OK key (Check mark)

Confirms the gray-highlighted menu items and scrolls forward in the menu.

Shows the first data record in the data log.

# Operation

## Switch on the instrument

Become familiar with the instrument functions before operation. Learn how to navigate through the menu and perform the relevant functions.

1. Connect the input connector of the power supply to the **POWER IN** port on the back of the TitrLab.
2. Connect the power supply to a power socket.
3. Press **ON/OFF** below the display to switch on the instrument.
4. The instrument performs a self-test upon each activation.

***Note:** Wait approximately 20 seconds before switching on again so as not to damage the electronics and mechanics of the instrument*

## Startup

### Language and data output setting

The titrator software supports multiple languages. The first time the instrument is switched on, a language selection list is displayed automatically upon completion of the self-test.

1. Select the required language.
2. Press **OK** to confirm the selection.
3. Select the type of data output (deactivated, printer, computer) and confirm with **OK**.

### Change the language setting

The instrument functions in the selected language until the option is changed.

1. Turn the instrument on.
2. Select **SYSTEM > LANGUAGE**.
3. Select the required language.
4. Press **OK** to confirm the selection.

### Change the data output setting

The instrument functions in the selected data output until the option is changed.

1. Turn the instrument on.
2. Select **SYSTEM > DATA OUTPUT**.
3. Select the required setting.
4. Press **OK** to confirm the setting.

## pH and Alkalinity content

1. Place a magnetic stirring bar in a titration cell and dispense the sample as follows:
  - Add 50 mL of the sample water.
2. Place the titration cell over the magnetic stirrer.
3. Lower the electrode into the titration cell.  
***Note:** Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure the combined pH Electrode, including the reference element porous pin is covered by sample.*
4. Select **ANALYZE PH + ALKALINITY**. "INSERT SAMPLE: 1" is shown on the display.
5. Select **START ANALYSIS**.  
The following steps are carried out automatically by TitraLab:
  - Measurement of pH
  - Subsequent commencement of titration
6. The result is shown on the display and, depending on the setting, sent to a printer or computer.
7. Press **OK** to confirm the result.
8. To begin another analysis, select **OK**. "NEXT SAMPLE" is shown on the display.

## pH measurement

1. Place a magnetic stirring bar and the sample in a titration cell.
2. Place the titration cell over the magnetic stirrer.
3. Lower the electrode into the titration cell.  
***Note:** Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure the combined pH Electrode, including the reference element porous pin is covered by sample.*
4. Select **MEASURE PH**. "STABILITY" or "CONTINUOUS" is shown on the display.
  - **STABILITY:** pH is shown on the display once stable.

- **CONTINUOUS:** pH is continuously shown on the display and, if a PC or printer is connected, transferred at a defined interval.
5. Select **START MEASURING**.
  6. The result is shown on the display and, depending on the setting, sent to a printer or computer.
  7. Press **OK** to confirm the result.
  8. To begin another analysis, select **OK**. "NEXT SAMPLE" is shown on the display.

## pH calibration

1. Place a magnetic stirring bar and pH 7.00 calibration solution in a titration cell.
2. Attach the titration cell to the magnetic stirrer.
3. Lower the electrode into the titration cell.  
***Note:** Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure the combined pH Electrode, including the reference element porous pin is covered by sample.*
4. Select **PH CALIBRATION**. "IMMERSE THE PH ELECTRODE IN THE FIRST BUFFER." is shown on the display.
5. Select **START**.
6. Following successful calibration, "NEXT POINT" or "FINISH CALIBRATION" is shown on the display.
  - **FINISH CALIBRATION:** Calibration finished as one-point calibration.
  - **NEXT POINT:** Place a magnetic stirring bar and pH 4.00 calibration solution in a titration cell. Continue as described.
7. "IMMERSE THE PH ELECTRODE IN THE SECOND BUFFER." is shown on the display.
8. Select **START**.
9. Following successful calibration, "FINISH CALIBRATION" is shown on the display.
10. Press **OK** to confirm.

## Manual component activation

### Burette

1. Select **MANUAL ACTIVATION > BURETTE** and press **OK** to confirm.
2. Select **FILL** with the **DOWN** navigation key to fill the burette.
3. Select **EMPTY** with the **UP** navigation key to empty the burette into the titration cell.

***Note:** Repeat this procedure (FILL/EMPTY) multiple times to make sure that the dispensing process functions correctly.*

4. Select **RETURN** with the **LEFT** navigation key to empty the burette back into the titrant bottle.

### Stirrer

1. Select **MANUAL ACTIVATION > STIRRER** and press **OK** to confirm.
2. Press **OK** to switch the stirrer on or off.
3. Use the navigation keys **UP** and **DOWN** to select the stirring speed.

## Maintenance

### WARNING

#### Risk of injury

Before maintenance or repair work is conducted, thoroughly rinse the instrument with deionized water and make sure that there are no residual chemicals in the tubes and pumps.

Only use genuine replacement parts and accessories.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

After maintenance tasks or service, the following points have to be checked:

- The housing of the instrument has to be closed correctly using the corresponding screws.

- Check the syringe cover for fissure or damage, replace if necessary. The syringe cover has to be attached firmly to the instrument using the corresponding screws.
- Place the instrument firmly on an even surface.

## Maintenance schedule

| Maintenance task                                          | Every week | 6 months | 12 months |
|-----------------------------------------------------------|------------|----------|-----------|
| Check all tube connections for leaks and damage.          | x          | x        |           |
| Replace burette syringe                                   |            | x        |           |
| Check valve block                                         |            |          | x         |
| Burette feed tube                                         |            |          | x         |
| Burette discharge tube                                    |            |          | x         |
| Bottle closure set                                        |            |          | x         |
| Check syringe and syringe connection for leaks and damage | x          |          |           |

# Troubleshooting

## Error messages

| Error displayed          | Cause                                                                                                      | Resolution                                                                                                                           |
|--------------------------|------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Maximum volume reached.  | The maximum volume is reached before the titration is completed. The titration is cancelled automatically. | Empty the titration cell and restart the titration.                                                                                  |
| Expired calibration.     | Validity of calibration has expired.                                                                       | Perform new calibration.                                                                                                             |
| Data logger empty.       | No results are stored or no results match the search filter.                                               | Change the search filter criteria.                                                                                                   |
| All data will be erased. | All data in the data log will be erased.                                                                   |                                                                                                                                      |
| Printing error.          | Printer is connected but not activated.                                                                    | Switch on the printer.                                                                                                               |
| PC not connected.        | Computer is activated but not connected.                                                                   | Check the computer connection.                                                                                                       |
| Unstable pH reading.     | Defective pH electrode and/or insufficient immersion.                                                      | Check electrode according to manufacturer recommendations, maintain electrode as required.<br>Check porous pin is covered by sample. |

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