DOC022.98.90360

TitraLab Salt analyzer

02/2013, Edition 2A

Основно ръководство за потребителя	bg
Základní uživatelská příručka	cs
Grundlæggende brugervejledning	da
Basis Bedienungsanleitung	de
Βασικό Εγχειρίδιο Χρήστη	el
Basic User Manual	en
Manual básico del usuario	es
Manuel d'utilisation de base	fr
Osnovni korisnički priručnik	hr
Alap felhasználói kézikönyv	hu
Manuale utente di base	it
Basisgebruikershandleiding	nl
Podstawowa instrukcja obsługi	pl
Manual do utilizador básico	pt
Manual de bază al utilizatorului	ro
Основное руководство пользователя	ru
Základný návod na použitie	sk
Osnovni uporabniški priročnik	sl
Grundläggande bruksanvisning	sv
Temel Kullanıcı Kılavuzu	tr

Български	3
Čeština	17
Dansk	31
Deutsch	45
Ελληνικά	59
English	75
Español	89
Français	103
Hrvatski	117
Magyar	131
Italiano	145
Nederlands	159
Polski	173
Português	187
Română	201
Русский	215
Slovenčina	229
Angleščina	243
Svenska	257
Türkçe	271

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Specifications

These are subject to change without notice!

Performance specifications		
Application area	LLV2303.97.1102	Determination of salt content in solid and liquid dairy products
	LLV2305.97.1102	Determination of salt content in tomato sauce
	LLV2307.97.1102	Determination of salt content in brine
	LLV2308.97.1102	Determination of salt content in water
Measuring range	LLV2303.97.1102	10 mL liquid samples 0,3–10 mL 1 g solid samples 0,25–5,6 % NaCl, 0,20–3,5 % Cl ⁻
	LLV2305.97.1102	1 g sample 0,25–5,6 % NaCl
	LLV2307.97.1102	0,5 mL sample 2,5–50 % NaCl
	LLV2308.97.1102	10 mL sample 10–350 ppm Cl ⁻
Resolution		Dependent on the selected units
Programs		3 titration programs
Languages		English, Spanish, French and Italian
Data logger		55 data records
Date/time		DD-MM-YYYY
Display		128 × 64 pt, backlit liquid crystal display
Keypad		7 membrane keys, splash-proof, PET, lifetime of >6 million key-presses
Burette	Resolution	0,001 mL
	Material of syringe	Polypropylene
	Volume	10 mL
Reagent tubes		White PTFE tubes
Peristaltic pump	Volumetric flow	70 mL/min
	Tube	Novoprene, internal diameter 2 mm)

Performance specifications		
Inputs and outputs	Measuring electrode	BNC connector
	Magnetic stirrer	RCA connector (cinch)
	Computer or printer	RJ11 telephone connector
	Power supply, external	24 VAC, 1,25 A, 30 VA
	Mouse	Mini-DIN connector (only available on LLV2307.97.1102 and LLV2308.97.1102)
Dimensions of measuring instrument and environmental conditions		
Enclosure material		ABS and enameled steel
Weight		Approximately 4 kg
Dimensions		130 × 160 × 300 mm
Working temperature		+15 to +40 °C (+60 to +104 °F)
Storage temperature		−10 to +50 °C (+14 to +122 °F)
Air humidity		<80 % relative humidity, non-condensing
Warranty		2 years

General information

Safety information

Please read the entire user manual carefully before the equipment is unpacked, set up or put into operation. Pay attention to all danger and warning statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

To make sure that the protection provided by this instrument is not impaired, do not use or install this instrument in any manner other than that specified in these operating instructions.

DANGER

Indicates a potentially or imminently hazardous situation that, if not avoided, results in death or serious injury.

WARNING

Indicates a potentially or imminently hazardous situation that, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation that, if not avoided, could result in damage to the instrument. Information that requires special emphasis.

Note: Information that supplements points in the main text.

Warning labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. For symbols applied to the instrument, corresponding warnings are provided in the user manual.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	This symbol may be found on an enclosure or barrier within the product and indicates a risk of electrical shock and/or death by electrocution.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems after August 12, 2005. In conformity with European local and national regulations, European electrical equipment users must now return old or end-of life equipment to the manufacturer for disposal at no charge to the user. Note: For return for recycling, please contact the equipment manufacturer or supplier for instructions on how to return end-of-life equipment, manufacturer-supplied electrical accessories, and all auxiliary items for correct disposal.

Chemical and biological safety

DANGER

Potential danger in the event of contact with chemical/biological materials.

Handling chemical samples, standards and reagents can be dangerous.

Familiarize yourself with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

Normal operation of this instrument may involve the use of hazardous chemicals or biologically harmful samples.

- Observe all cautionary information printed on the original solution containers and safety data sheet prior to their use.
- Dispose of all consumed solutions in accordance with national regulations and laws.
- Select the type of protective equipment suitable to the concentration and quantity of the dangerous material at the respective work place.

Installation

WARNING

Risk of injury to eyes, skin and respiratory

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Make sure that the system is only unpacked, assembled, connected and operated by qualified and trained personnel.


Unpack the instrument

Unpack all supplied parts carefully, as they are highly sensitive in part to shock and impact. Read the user manual prior to installation and proceed exactly as described.

Scope of delivery

Please check that the order is complete. If anything is missing or damaged, please contact the manufacturer or retailer.

Figure 1 Scope of delivery

1	TitraLab salt system	8	Silver electrode with cable and tube guide
2	2× tube with bottle closure (DIN 45 thread) supplied pre-installed	9	Mouse (versions LLV2307.97.1102 and LLV2308.97.1102 only)
3	2× tube to attach to electrode, supplied pre-installed	10	Adapter for power supply
4	Titration stand holder including magnetic stirrer	11	Power supply, external
5	10×white O-rings for burette and pump threads	12	Hardware sleeve
6	10× Teflon magnetic stirring bars (20 mm × 6 mm)	13	Syringe with plunger
7	10×titration cells, 50 mL, Ø 30 mm, polypropylene with lid	14	Desiccant cartridge

Operating environment

Note the following points so that the instrument will operate perfectly and have a long service life.

- Place the instrument firmly on an even surface. Do not push any objects under the instrument.
- The ambient temperature must be 15 to 40 °C (60 to 104 °F).
- The relative humidity should be less than 80 %; moisture should not condense on the instrument.
- A minimum clearance of 15 cm (5,9 in.) must be ensured above and on all sides of the device; this allows the air to circulate and prevents overheating of the electric parts. The dosage arm of the stepper motor operating the syringe piston must be able to move freely.
- Do not operate or store the instrument in extremely dusty, humid or wet locations.
- Make sure that no liquids enter the instrument and immediately wipe away any liquids that contact the instrument.
- Protect the instrument against vibrations, direct sunlight, corrosive gases as well as strong magnetic and/or electric fields.
- Only use genuine replacement parts and accessories.
- The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

Installation

Installation of the burette

The burette dispenses the titrant, a Silver Nitrate (AgNO_3) solution that, depending on the area of application, is either 0,01 M, 0,1 M or 0,5 M.

The burette comprises the syringe, stepper motor and valve. The stepper motor and valve are already installed in the TitraLab; the syringe has to be installed by the user.

WARNING

Risk of injury to pinch points and cuts

Make sure that the titrator is never used/operated without the syringe cover.

Make sure that the system is only assembled by qualified and trained personnel.

NOTICE

Warning against system damage

Make sure that the syringe is never operated without the attachment sleeve.

Only use genuine replacement parts and accessories.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

Note: *Make sure that you do not misplace the attachment sleeve for the syringe. It is not provided with a new syringe.*

1. Push the syringe piston into the syringe body.
Slide the attachment sleeve fully onto the syringe so that the narrow end points toward the piston.
2. Remove the syringe cover by unscrewing it.
3. Point the syringe tip upward, place the syringe into the bracket such that the attachment sleeve is above the bracket.
4. With the narrow side facing down, carefully press the attachment sleeve into the bracket.
Hold the attachment sleeve in the bracket and screw the syringe onto the valve from below.
5. Withdraw the piston from the syringe and place it onto the screw in the dosage arm.
Hold the syringe piston and turn the attachment screw on the dosage arm into the thread in the piston.
6. Fix the syringe cover with six screws.


⚠ WARNING

Risk of injury to eyes, skin and respiratory

During the first days of system use, check the screw joint between the syringe and valve to avoid leaks and the potential escape of reagents.

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Installation of the burette

Installation of the titration stand

The titration stand consists of a holder onto which the electrode, titration cell and magnetic stirrer are attached. The titration is performed in the titration cell.

1. Place the titration stand on the right side (syringe side) of the titrator (observe the routing of the cable).
2. Place a magnetic stirring bar into a titration cell. Attach the titration cell to the magnetic stirrer.
3. Place the rubber ring on the stand rod to set the minimum height of the electrode holder.
4. Press the button to lock the electrode holder and place the electrode holder on the stand rod. The electrode holder is locked upon release of the button.

5. Remove the white cap from the tip of the electrode. Insert the electrode into the electrode holder.

NOTICE


Warning against system damage

Make sure that the height setting of the electrode holder protects the electrode against damage from the stirring bar and against contact with the bottom of the titration cell.

Keep the lock button pressed whenever the electrode holder is being moved.

Note: Whenever the electrode is taken from the titration stand or the titration cell is removed from the titration stand, always place a container beneath the electrode. Otherwise, drips from the sample or chemicals could contaminate the work area.

Installation of the titration stand

Electrical connection of the magnetic stirrer

1. Connect the cable of the magnetic stirrer to the **STIRRER** port on the back of the TitraLab.
The power is supplied via the TitraLab.

Installation of tubes

WARNING

Risk of injury to eyes, skin and respiratory

During the first days of system use, check the screw joint between the syringe and valve to avoid leaks and the potential escape of reagents.

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Make sure that the tubes are never kinked.

Regularly check the system as well as all tubes and connections for leaks, condition and cleanliness.


To prevent twisting of the tubes, you should always connect flange tube connections before those with grommets.

The system is delivered with four tubes; these are supplied pre-installed on the TitraLab and only need to be connected to the relevant reagent bottles and to the electrode.

Removal of protective stickers (Uninstalled tubes configuration)

1. Remove the protective stickers on the burette valve and the four pump connections.
2. Make sure that the pre-installed white O-rings are not removed.
3. Check that the O-rings are in a horizontal position inside the thread.

Figure 2 Position of protective stickers


Connection of feed tubes

The system is delivered with two feed tubes with bottle closures (DIN 45 thread). These feed tubes are supplied pre-installed on the TitraLab and only need to be connected to the relevant reagent bottles.

1. Connect the tube from the IN socket of the burette valve to the titrant bottle (AgNO_3). Screw on the bottle closure to do this.
2. Connect the tube from the IN socket of the pump to the acid bottle (HNO_3). Screw on the bottle closure to do this.
3. Place a desiccant cartridge into the designated opening in each bottle closure.

Figure 3 Connection of feed tubes

Connection of the electrode

Electrical connection of the electrode


1. Connect the electrode cable to the **ELECTRODE** port on the back of the TitraLab.
The power is supplied via the TitraLab.

Connection of the electrode tubes

The sensor holder of the titration stand has two tube guides.

1. Connect the tube from the OUT socket of the pump to the electrode.
2. To do this, push the free end of the tube into guide "1" of the sensor holder until it reaches the middle of the electrode.
3. Connect the tube on the OUT socket of the burette valve to the electrode.
4. To do this, push the free end of the tube into guide "2" of the sensor holder until the end of the tube is flush with the tip of the electrode.

Figure 4 Connection of the electrode

Power supply, external

⚠ DANGER


Risk of electric shock

Make sure the AC supply voltage is suitable for the power supply and the mains plug style is suitable for the type of outlet socket.

Only use the original power supply and plug connectors.

The power supply is provided with different plug connectors. The correct plug connector must be attached to the power supply before first use.

Figure 5 Plug connector attachment

Note: To change the plug connector, press in the catch and slide the plug connector upward.

Keypad

Navigation key UP

- Used to select between different options on the display.
- Increases numeric values.
- Shows the previous data record in the data log menu.
- Empties the syringe through the discharge tube during manual open loop control.

Navigation key DOWN

- Used to select between different options on the display.
- Reduces numeric values.
- Shows the next data record in the data log menu.
- Fills the syringe through the feed tube during manual open loop control.

ON/OFF key

Navigation key LEFT

- Moves the cursor backward for entry of numeric values.
- Scrolls up in menus.
- Empties the syringe through the feed tube during manual open loop control.

Navigation key RIGHT

- Moves the cursor forward for entry of numeric values.
- Scrolls down in menus.

RETURN key

- Scrolls back in menus.
- Cancels a manual control.
- Cancels a titration.

OK key (Check mark)

- Confirms the gray-highlighted menu items and scrolls forward in the menu.
- Shows the first data record in the data log.

Operation

Switch on the instrument

Become familiar with the instrument functions before operation. Learn how to navigate through the menu and perform the relevant functions.

1. Connect the input connector of the power supply to the **POWER IN** port on the back of the TitraLab.
2. Connect the power supply to a power socket.
3. Press **ON/OFF** below the display to switch on the instrument.
4. The instrument performs a self-test upon each activation.

Note: Wait approximately 20 seconds before switching on again so as not to damage the electronics and mechanics of the instrument

Startup

Language and data output setting

The titrator software supports multiple languages. The first time the instrument is switched on, a language selection list is displayed automatically upon completion of the self-test.

1. Select the required language.
2. Press **OK** to confirm the selection.
3. Select the type of data output (deactivated, printer, computer) and confirm with **OK**.

Change the language setting

The instrument functions in the selected language until the option is changed.

1. Turn the instrument on.
2. Select **SYSTEM > LANGUAGE**.
3. Select the required language.
4. Press **OK** to confirm the selection.

Chloride analysis in solid samples

1. Place a magnetic stirring bar in a titration cell and dispense the sample as follows:
 - **Cheese or solid dairy products:** Add 1 g of the sample dissolved in approximately 10 mL lukewarm distilled water
 - **Tomato sauces:** Add 1 g of the sample dissolved in approximately 10 mL of distilled water
 - **Brine:** Add 0,5 g of the sample dissolved in approximately 10 mL of distilled water
2. Attach the titration cell to the magnetic stirrer.
3. Before the analysis, stir the sample for at least 30 seconds in order to extract the salt in full.
4. Lower the electrode into the titration cell.

Note: Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure the combined Silver rod electrode, including the reference element porous pin is covered by sample.
5. Select **ANALYZE CHLORIDE (W)**. "**SAMPLE: 1 WEIGHT**" is shown on the display.
6. Enter the weight of the sample and confirm with **OK**.

The following steps are carried out automatically by TitraLab:

 - Acid is dispensed into the titration cell for 2 seconds
 - The solution is stirred for 3 seconds
 - Titration is then started
7. The result is shown on the display and, depending on the setting, sent to a printer or computer.
8. Press **OK** to confirm the result.
9. To begin another analysis, select **OK**. "**NEXT SAMPLE**" is shown on the display.

Chloride analysis in liquid samples

1. Place a magnetic stirring bar in a titration cell and dispense the sample as follows:
 - **Milk or liquid dairy products:** 10 mL sample
 - **Water:** 10 mL sample
 2. Attach the titration cell to the magnetic stirrer.
 3. Lower the electrode into the titration cell.

Note: Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure the combined Silver rod electrode, including the reference element porous pin is covered by sample.
 4. Select **ANALYZE CHLORIDE (V)**. "INSERT SAMPLE: 1" is shown on the display.
 5. Select **OK**.

The following steps are carried out automatically by TitraLab:

 - Acid is dispensed into the titration cell for 2 seconds
 - The solution is stirred for 3 seconds
 - Titration is then started
 6. The result is shown on the display and, depending on the setting, sent to a printer or computer.
 7. Press **OK** to confirm the result.
- To begin another analysis, select **OK**. "NEXT SAMPLE" is shown on the display.

Manual component activation

Burette

1. Select **MANUAL ACTIVATION > BURETTE** and press **OK** to confirm.
 - Select **FILL** with the **DOWN** navigation key to fill the burette.
 - Select **EMPTY** with the **UP** navigation key to empty the burette into the titration cell.
 - Select **RETURN** with the **LEFT** navigation key to empty the burette back into the titrant bottle.

Note: Repeat this procedure multiple times to make sure that the dispensing process functions correctly.

Pump

1. Select **MANUAL ACTIVATION > PUMPS** and press **OK** to confirm.
2. Press **OK** to switch the pump on or off.

Note: Repeat this procedure multiple times to make sure that the dispensing process functions correctly.

Stirrer

1. Select **MANUAL ACTIVATION > STIRRER** and press **OK** to confirm.
2. Press **OK** to switch the stirrer on or off.
3. Use the navigation keys **UP** and **DOWN** to select the stirring speed.

Maintenance

WARNING

Risk of injury

Before maintenance or repair work is conducted, thoroughly rinse the instrument with deionized water and make sure that there are no residual chemicals in the tubes and pumps.

Only use genuine replacement parts and accessories.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

Maintenance schedule

Maintenance task	Every week	6 months	12 months
Replace dry desiccant cartridges content.	x		
Check all tube connections for leaks and damage.	x		
Replace burette syringe		x	
Valve block			x
Burette feed tube			x
Burette discharge tube			x
Peristaltic pump, pump cassette			x
Pump feed tube			x
Pump discharge tube			x
Bottle closure set			x

Troubleshooting

Error messages

Error displayed	Cause	Resolution
Maximum volume reached.	The maximum volume is reached before the titration is completed. The titration is cancelled automatically.	Empty the titration cell and restart the titration.
Equivalence points not found.	No equivalence points found during titration.	Check system with Sodium chloride standard of known concentration. Check electrode according to manufacturer recommendations, maintain electrode as required. Check porous pin is covered by sample.
	Chloride content below measuring range.	
Data logger empty.	No results are stored or no results match the search filter.	Change the search filter criteria.
All data will be erased.	All data in the data log will be erased.	
Printing error.	Printer is connected but not activated.	Switch on the printer.
PC not connected.	Computer is activated but not connected.	Check the computer connection.
Unstable mV reading	Defective Silver rod electrode and/or insufficient immersion.	Check electrode according to manufacturer recommendations, maintain electrode as required. Check porous pin is covered by sample.

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80539-0389 U.S.A.
(U.S.A. only)