Sonda 50 44

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

EN

#### 50 44 probe

Refer to the meter manual for more information.

#### Caution

⚠ Personal injury hazard. Glass components can break. Handle with care to prevent cuts.

⚠ Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.

#### Introduction

This refillable reference electrode uses a lithium acetate electrolyte.

#### Technical information

|                |                                   |
|----------------|-----------------------------------|
| Junction       | Ceramic                           |
| Reference      | Ag/AgCl and salt bridge           |
| Operating temp | 0 to 80 °C (32 to 176 °F)         |
| Limits         | Colloids or solids in suspension. |

#### Notice

Always use the reference electrode with an indicator electrode.

#### 1 Connect the probe

#### 2 Preparation

#### 3 Calibration/Measurement

Rinse the probe with deionized water and store in solution between measurements.

#### 4 Maintenance

Refill the probe as needed to keep it full. Keep the junction clean. Replace if damaged or dry. After replacement, refill the electrode and soak junction in LiCH<sub>3</sub>COO 0.1M solution for 30 minutes.

Cleaning: soak the probe in a cleaning solution. Rinse with deionized water.

Hydration: soak the dry probe in a dilute HCl solution for several hours.

#### 5 Storage

Put the probe in the protector with LiCH<sub>3</sub>COO 0.1M solution.

#### Troubleshooting

Be sure that the samples are suitable for the application and that the junction is clean.

#### Warranty

1 year for manufacturing faults only. Damage from use is not covered.

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