DOC022.98.90359

TitraLab SO₂ analyzer

02/2013, Edition 2A

Основно ръководство за потребителя	bg
Základní uživatelská příručka	cs
Grundlæggende brugervejledning	da
Basis Bedienungsanleitung	de
Βασικό Εγχειρίδιο Χρήστη	el
Basic User Manual	en
Manual básico del usuario	es
Manuel d'utilisation de base	fr
Osnovni korisnički priručnik	hr
Alap felhasználói kézikönyv	hu
Manuale utente di base	it
Basisgebruikershandleiding	nl
Podstawowa instrukcja obsługi	pl
Manual do utilizador básico	pt
Manual de bază al utilizatorului	ro
Основное руководство пользователя	ru
Základný návod na použitie	sk
Osnovni uporabniški priročnik	sl
Grundläggande bruksanvisning	sv
Temel Kullanıcı Kılavuzu	tr

Български	3
Čeština	19
Dansk	35
Deutsch	49
Ελληνικά	63
English	79
Español	93
Français	109
Hrvatski	125
Magyar	141
Italiano	157
Nederlands	173
Polski	189
Português	205
Română	221
Русский	237
Slovenčina	253
Angleščina	269
Svenska	283
Türkçe	297

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Specifications

These are subject to change without notice!

Performance specifications		
Measuring principle		Potentiometric titration according to the Ripper method for the determination of free SO ₂ and total SO ₂ content in wine
Measuring range		0–640 mg/L SO ₂
Resolution		1 mg/L SO ₂
Programs		3 titration programs
Languages		English, Spanish, French and Italian
Data logger		55 data records
Date/time		DD-MM-YYYY
Display		128 × 64 pt, backlit liquid crystal display
Keypad		7 membrane keys, splash-proof, PET, lifetime of >6 million key-presses
Burette	Resolution	0,001 mL
	Material of syringe	Polypropylene
	Volume	10 mL
Reagent tubes		White PTFE tubes
Peristaltic pump	Volumetric flow	70 mL/min
	Tube	Novoprene, internal diameter 2 mm

Performance specifications		
Inputs and outputs	Measuring electrode	BNC connector
	Magnetic stirrer	RCA connector (cinch)
	Computer or printer	RJ11 telephone connector
	Power supply, external	24 VAC, 1,25 A, 30 VA
	Mouse	Mini-DIN connector
Dimensions of measuring instrument and environmental conditions		
Enclosure material		ABS and enameled steel
Weight		approximately 4 kg
Dimensions		130 × 160 × 300 mm
Working temperature		+15 to +40 °C (+60 to +104 °F)
Storage temperature		-10 to +50 °C (+14 to +122 °F)
Air humidity		<80 % relative humidity, non-condensing
Warranty		
Warranty		2 years

General information

Safety information

Please read the entire user manual carefully before the equipment is unpacked, set up or put into operation. Pay attention to all danger and warning statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

To make sure that the protection provided by this instrument is not impaired, do not use or install this instrument in any manner other than that specified in these operating instructions.

DANGER

Indicates a potentially or imminently hazardous situation that, if not avoided, results in death or serious injury.

WARNING

Indicates a potentially or imminently hazardous situation that, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation that, if not avoided, could result in damage to the instrument. Information that requires special emphasis.

Note: Information that supplements points in the main text.

Warning labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. For symbols applied to the instrument, corresponding warnings are provided in the user manual.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	This symbol may be found on an enclosure or barrier within the product and indicates a risk of electrical shock and/or death by electrocution.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems after August 12, 2005. In conformity with European local and national regulations, European electrical equipment users must now return old or end-of life equipment to the manufacturer for disposal at no charge to the user. Note: For return for recycling, please contact the equipment manufacturer or supplier for instructions on how to return end-of-life equipment, manufacturer-supplied electrical accessories, and all auxiliary items for correct disposal.

Chemical and biological safety

DANGER

Potential danger in the event of contact with chemical/biological materials.

Handling chemical samples, standards and reagents can be dangerous.

Familiarize yourself with the necessary safety procedures and the correct handling of the chemicals before the work and read and follow all relevant safety data sheets.

Normal operation of this instrument may involve the use of hazardous chemicals or biologically harmful samples.

- Observe all cautionary information printed on the original solution containers and safety data sheet prior to their use.
- Dispose of all consumed solutions in accordance with national regulations and laws.
- Select the type of protective equipment suitable to the concentration and quantity of the dangerous material at the respective work place.

Installation

WARNING

Risk of injury to eyes, skin and respiratory

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Make sure that the system is only unpacked, assembled, connected and operated by qualified and trained personnel.


Unpack the instrument

Unpack all supplied parts carefully, as they are highly sensitive in part to shock and impact. Read the user manual prior to installation and proceed exactly as described.

Scope of delivery

Please check that the order is complete. If anything is missing or damaged, please contact the manufacturer or retailer.

Figure 1 Scope of delivery

1	TitraLab SO ₂ titrator	8	Mouse
2	3× tube with bottle closure (DIN 45 thread) supplied pre-installed	9	Power supply, external
3	Titration stand with magnetic stirrer	10	Adapter for power supply
4	10× titration cells, 50 mL, Ø 30 mm, polypropylene with lid	11	Hardware sleeve
5	10× Teflon magnetic stirring bars (20 mm × 6 mm)	12	Syringe with plunger
6	10× white O-rings for burette and pump threads	13	3× desiccant cartridge
7	Double-platinum electrode including cable and tubes		

Operating environment

Note the following points so that the instrument will operate correctly and have a long service life.

- Place the instrument firmly on an even surface. Do not push any objects under the instrument.
- The ambient temperature must be 15 to 40 °C (60 to 104 °F).
- The relative humidity should be less than 80 %; moisture should not condense on the instrument.
- A minimum clearance of 15 cm (5,9 in.) must be ensured above and on all sides of the device; this allows the air to circulate and prevents overheating of the electric parts. The dosage arm of the stepper motor operating the syringe piston must be able to move freely.
- Do not operate or store the instrument in extremely dusty, humid or wet locations.
- Make sure that no liquids enter the instrument and immediately wipe away any liquids that contact the instrument.
- Protect the instrument against vibrations, direct sunlight, corrosive gases as well as strong magnetic and/or electric fields.
- Only use genuine replacement parts and accessories.
- The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

Installation

Installation of the burette

The burette dispenses the titrant, i.e. the 0,01 M I_3 solution.

The burette comprises the syringe, stepper motor and valve. The stepper motor and valve are already installed in the TitraLab; the syringe must be installed by the user.

WARNING

Risk of injury to pinch points and cuts

Make sure that the titrator is never used/operated without the syringe cover.

Make sure that the system is only assembled by qualified and trained personnel.

NOTICE

Warning against system damage

Make sure that the syringe is never operated without the attachment sleeve.

Only use genuine replacement parts and accessories.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

Note: Make sure that you do not misplace the attachment sleeve for the syringe. It is not provided with a new syringe.

1. Push the syringe piston into the syringe body.
Slide the attachment sleeve fully onto the syringe so that the narrow end points toward the piston.
2. Remove the syringe cover by unscrewing it.
3. Point the syringe tip upward, place the syringe into the bracket such that the attachment sleeve is above the bracket.
4. With the narrow side facing down, carefully press the attachment sleeve into the bracket.
Hold the attachment sleeve in the bracket and screw the syringe onto the valve from below.
5. Withdraw the piston from the syringe and place it onto the screw in the dosage arm.
Hold the syringe piston and turn the attachment screw on the dosage arm into the thread in the piston.
6. Fix the syringe cover with six screws.


⚠ WARNING

Risk of injury to eyes, skin and respiratory

During the first days of system use, check the screw joint between the syringe and valve to avoid leaks and the potential escape of reagents.

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Figure 2 Installation of the burette

Installation of the titration stand

The titration stand consists of a holder onto which the electrode, titration cell and magnetic stirrer are attached. The titration is performed in the titration cell.

1. Place the titration stand on the right side (syringe side) of the titrator (observe the routing of the cable).
2. Place a magnetic stirring bar into a titration cell. Attach the titration cell to the magnetic stirrer.
3. Place the rubber ring on the stand rod to set the minimum height of the electrode holder.
4. Press the button to lock the electrode holder and place the electrode holder on the stand rod. The electrode holder is locked upon release of the button.

5. Remove the white cap from the tip of the electrode. Insert the electrode into the electrode holder.

NOTICE


Warning against system damage

Make sure that the height setting of the electrode holder protects the electrode against damage from the stirring bar and against contact with the bottom of the titration cell.

Keep the lock button pressed whenever the electrode holder is being moved.

Note: Whenever the electrode is taken from the titration stand or the titration cell is removed from the titration stand, always place a container beneath the electrode. Otherwise, drips from the sample or chemicals could contaminate the work area.

Figure 3 Installation of the titration stand

Electrical connection of the magnetic stirrer

1. Connect the cable of the magnetic stirrer to the **STIRRER** port on the back of the TitrLab.
The power is supplied via the TitrLab.

Installation of tubes

WARNING

Risk of injury to eyes, skin and respiratory

During the first days of system use, check the screw joint between the syringe and valve to avoid leaks and the potential escape of reagents.

When working with chemicals and/or solvents, always observe the relevant accident prevention regulations and wear appropriate personal protective clothing for eye, face, hand, body and/or respiratory protection.

Make sure that the tubes are never kinked.


Regularly check the system as well as all tubes and connections for leaks, condition and cleanliness.

To prevent twisting of the tubes, you should always connect flange tube connections before those with grommets.

Removal of protective stickers (Uninstalled tubes configuration)

1. Remove the protective stickers on the burette valve and the four pump connections.
2. Make sure that the pre-installed white O-rings are not removed.
3. Check that the O-rings are in a horizontal position inside the thread.

Figure 4 Position of protective stickers


Connection of feed tubes

The system is delivered with three feed tubes with bottle closures (DIN 45 thread). These feed tubes are supplied pre-installed on the TitrLab and only need to be connected to the relevant reagent bottles.

1. Connect the tube from the burette valve to the titrant bottle. Screw onto the bottle closure to do this.
2. Connect the tube from the IN socket of the **ACID** pump to the acid bottle. Screw on the bottle closure to do this.
3. Connect the tube from the IN socket of the **ALKALI** pump to the alkali bottle. Screw on the bottle closure to do this.
4. Place a desiccant cartridge in each bottle closure.

Figure 5 Connection of feed tubes

Connection of the electrode


1. Connect the electrode cable to the **SO₂-ELECTRODE** port on the back of the TitraLab.
The power is supplied via the TitraLab.

Connection of the electrode tubes

The electrode has three tubes with screw connections.

1. Screw one tube onto the burette valve.
2. Screw one tube into the OUT socket of the **ACID** pump.
3. Screw one tube into the OUT socket on the **ALKALI** pump.

Figure 6 Connection of the electrode

Power supply, external

⚠ DANGER


Risk of electric shock

Make sure the AC supply voltage is suitable for the power supply and the mains plug style is suitable for the type of outlet socket.

Only use the original power supply and plug connectors.

The power supply is provided with different plug connectors. The correct plug connector must be attached to the power supply before first use.

Figure 7 Plug connector attachment

Note: To change the plug connector, press in the catch and slide the plug connector upward.

Keypad

Navigation key UP

- Used to select between different options on the display.
- Increases numeric values.
- Shows the previous data record in the data log menu.
- Empties the syringe through the discharge tube during manual open loop control.

Navigation key DOWN

- Used to select between different options on the display.
- Reduces numeric values.
- Shows the next data record in the data log menu.
- Fills the syringe through the feed tube during manual open loop control.

ON/OFF key

Navigation key LEFT

- Moves the cursor backward for entry of numeric values.
- Scrolls up in menus.
- Empties the syringe through the feed tube during manual open loop control.

Navigation key RIGHT

- Moves the cursor forward for entry of numeric values.
- Scrolls down in menus.

RETURN key

- Scrolls back in menus.
- Cancels a manual control.
- Cancels a titration.

OK key (Check mark)

- Confirms the gray-highlighted menu items and scrolls forward in the menu.
- Shows the first data record in the data log.

Operation

Switch on the instrument

Become familiar with the instrument functions before operation. Learn how to navigate through the menu and perform the relevant functions.

1. Connect the input connector of the power supply to the **POWER IN** port on the back of the TitraLab.
2. Connect the power supply to a power socket.
3. Press **ON/OFF** below the display to switch on the instrument.
4. The instrument performs a self-test upon each activation.

Note: Wait approximately 20 seconds before switching on again so as not to damage the electronics and mechanics of the instrument

Startup

Language selection and data output setting

The titrator software supports multiple languages. The first time the instrument is switched on, a language selection list is displayed automatically upon completion of the self-test.

1. Select the required language.
2. Press **OK** to confirm the selection.
3. Select the type of data output (deactivated, printer, computer) and confirm with **OK**.

Change the language setting

The instrument functions in the selected language until the option is changed.

1. Turn the instrument on.
2. Select **SYSTEM > LANGUAGE**.
3. Select the required language.
4. Press **OK** to confirm the selection.

Analysis of free SO₂

1. Dispense 10 mL of the sample into a titration cell.

Note: If the sample is shaken or filtered, the readings will be lower.

Note: Work with a sample and room temperature of (20 °C [70 °F]). Higher temperatures lead to higher readings; lower temperatures lead to lower readings.

2. Obtain just enough solid potassium iodide (KI) to cover the tip of a spoon (approximately 0,2 g) and add to the titration cell.
3. Add a magnetic stirring bar and attach the titration cell to the magnetic stirrer.
4. Lower the electrode into the titration cell until the **"FREE"** mark on the electrode is flush with the top edge of the titration cell.

Note: Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure that the electrode wires are covered by sample.

5. Select **ANALYZE FREE SO₂**. **"INSERT SAMPLE: 1"** is shown on the display.
6. Select **START ANALYSIS**.

The following steps are carried out automatically by TitraLab:

- Acid is dispensed into the titration cell for 1 second
 - Titration is then started
7. The result is shown on the display and, depending on the setting, sent to a printer or computer.
 8. Press **OK** to confirm the result.
 9. To begin another analysis, select **OK**. **"NEXT SAMPLE"** is shown on the display.

Analysis of total SO₂

1. Dispense 10 mL of the sample into a titration cell.
Note: *If the sample is shaken or filtered, the readings will be lower.*
Note: *Work with a sample and room temperature of (20 °C [70 °F]). Higher temperatures lead to higher readings; lower temperatures lead to lower readings.*
2. Obtain just enough solid potassium iodide (KI) to cover the tip of a spoon (approximately 0,2 g) and add to the titration cell.
3. Add a magnetic stirring bar and attach the titration cell to the magnetic stirrer.
4. Lower the electrode into the titration cell until the **"TOTAL"** mark on the electrode is flush with the top edge of the titration cell.
Note: *Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure that the electrode wires are covered by sample.*
5. Select **ANALYZE TOTAL SO2**. **"INSERT SAMPLE: 1"** is shown on the display.
6. Select **START ANALYSIS**.
The following steps are carried out automatically by TitraLab:
 - Alkali is dispensed into the titration cell for 3 seconds
 - The solution is stirred for 3 seconds
 - The reaction time is 5 minutes; a countdown on the display indicates the time remaining
 - Acid is dispensed into the titration cell for 5 seconds
 - Titration is then started
7. The result is shown on the display and, depending on the setting, sent to a printer or computer.
8. Press **OK** to confirm the result.
9. To begin another analysis, select **OK**. **"NEXT SAMPLE"** is shown on the display.

Analysis of total SO₂ in a series of samples

A reaction time of 5 minutes is required to determine the total SO₂ content. TitraLab is locked and unavailable for other analyses during this time. If the total SO₂ is to be determined in multiple samples, the best program to select is **"PREPARE + ANALYZE"**.

1. Dispense 10 mL of each sample into each titration cell of your series.
Note: *If the sample is shaken or filtered, the readings will be lower.*
Note: *Work with a sample and room temperature of (20 °C [70 °F]). Higher temperatures lead to higher readings; lower temperatures lead to lower readings.*
2. Obtain just enough solid potassium iodide (KI) to cover the tip of a spoon (approximately 0,2 g) and add to each titration cells.
3. Add a magnetic stirring bar to each one.
4. Attach the titration cell with the first sample to the magnetic stirrer.
5. Lower the electrode into the titration cell until the **"TOTAL"** mark on the electrode is flush with the top edge of the titration cell.
Note: *Make sure that the electrode cannot be damaged by the stirring bar or come into contact with the bottom of the titration cell. Make sure that the electrode wires are covered by sample.*
6. Select **PREPARE + ANALYZE**. **"INSERT SAMPLE: 1"** is shown on the display.
7. Select **START**.
The following steps are carried out automatically by TitraLab:
 - Alkali is dispensed into the titration cell for 3 seconds
 - The solution is stirred for 3 seconds
 - The reaction time is 5 minutes; a countdown on the display indicates the time remaining
8. Remove the titration cell with the first sample and attach the titration cell with the second sample to the magnetic stirrer.
9. Lower the electrode into the titration cell until the **"TOTAL"** mark on the electrode is flush with the top edge of the titration cell.

10. Select **NEXT SAMPLE**.

The following steps are carried out automatically by TitraLab:

- Alkali is dispensed into the titration cell for 3 seconds
- The solution is stirred for 3 seconds

11. Process all other samples in the series in the same way.

12. Five signal tones indicate the progression of the reaction time (5 minutes) for the 1st sample in the series.

Decide whether to stop preparing samples and start running analyses of prepared samples or continue preparing sequence.

13. If not all the samples in the series have been prepared, select **NEXT SAMPLE**.

14. Once all samples have been prepared, attach the titration cell with the first sample to the magnetic stirrer.

15. Lower the electrode into the titration cell until the "**TOTAL**" mark on the electrode is flush with the top edge of the titration cell.

16. Select **START ANALYSIS**.

The following steps are carried out automatically by TitraLab:

- Acid is dispensed into the titration cell for 5 seconds
- Titration is then started

17. The result is shown on the display and, depending on the setting, sent to a printer or computer.

18. Press **OK** to confirm the result.

19. Remove the titration cell with the first sample and attach the titration cell with the second sample to the magnetic stirrer.

20. Lower the electrode into the titration cell until the "**TOTAL**" mark on the electrode is flush with the top edge of the titration cell.

21. Select **NEXT SAMPLE**. The second sample is titrated. Process all subsequent samples in the series in the same way.

Manual component activation

Burette

1. Select **MANUAL ACTIVATION > BURETTE** and press **OK** to confirm.
 - Select **FILL** with the **DOWN** navigation key to fill the burette.
 - Select **EMPTY** with the **UP** navigation key to empty the burette into the titration cell.
 - Select **RETURN** with the **LEFT** navigation key to empty the burette back into the titrant bottle.

Note: Repeat this procedure multiple times to make sure that the dispensing process functions correctly.

Pumps

1. Select **MANUAL ACTIVATION > PUMPS ACTIVATION**, and press **OK** to confirm.
2. **PUMP1 DISPENSING** selects the acid pump.
PUMP2 DISPENSING selects the alkali pump.
CONTINUE... activates the selected pump.
3. Press **OK** to switch the pump on or off.

Note: Repeat this procedure multiple times to make sure that the dispensing process functions correctly.

Stirrer

1. Select **MANUAL ACTIVATION > STIRRER** and press **OK** to confirm.
2. Press **OK** to switch the stirrer on or off.
3. Use the navigation keys **UP** and **DOWN** to select the stirring speed.

Maintenance

WARNING

Risk of injury

Before maintenance or repair work is conducted, thoroughly rinse the instrument with deionized water and make sure that there are no residual chemicals in the tubes and pumps.

Only use genuine replacement parts and accessories.

The instrument must be serviced by the manufacturer's service department at the prescribed intervals.

Maintenance schedule

Maintenance task	Every week	6 months	12 months
Replace dry desiccant cartridges content	x		
Check all tube connections for leaks and damage	x		
Replace burette syringe		x	
Valve block			x
Burette feed tube			x
Burette discharge tube			x
Peristaltic pump, pump cassette			x
Pump feed tube			x
Pump discharge tube			x
Bottle closure set			x

Troubleshooting

Error messages

Error displayed	Cause	Resolution
Maximum volume reached.	The maximum volume is reached before the titration is completed. The titration is cancelled automatically.	Empty the titration cell and restart the titration.
Error initial tendency.	The first measurement is incorrect.	When measuring free SO ₂ : The sample might not contain any SO ₂ . When measuring total SO ₂ : Make sure that the acid and alkali pumps are functional.
Data logger empty.	No results are stored or no results match the search filter.	Change the search filter criteria.
All data will be erased.	All data in the data log will be erased.	
Printing error.	Printer is connected but not activated.	Switch on the printer.
PC not connected.	Computer is activated but not connected.	Check the computer connection.

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80539-0389 U.S.A.
(U.S.A. only)