

52 11

DOC022.98.90540
06/2016, Edition 2

Puntos clave. *Punti chiave.* Points clé. *Key points.*

- ① Tapón protector.
Tappo protettore.
Bouchon protecteur.
Protective cup.
- ② Cabezal con conector S7
Testa con connettore S7
Tête avec connecteur S7
S7 connector head
- ③ Orificio de relleno de electrolito.
Orificio di riempimento per l'elettrolita.
Orifice de remplissage de l'électrolyte.
Electrolyte filling port.
- ④ Electrolito de referencia.
Elettrolita di riferimento.
Électrolyte de référence.
Reference electrolyte.
- ⑤ Protector de almacenamiento.
Protettore per la conservazione.
Protecteur de stockage.
Storage protector.
- ⑥ Elemento de referencia.
Elemento di riferimento.
Élément de référence.
Reference element.
- ⑦ Diafragmas cerámicos (3).
Diaframmi ceramici (3).
Diaphragmes en céramique (3).
Ceramic Diaphragms (3).
- ⑧ Membrana sensible.
Membrana sensibile.
Membrane sensible.
Sensitive glass membrane.

Electrodo de pH 52 11

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pH electrode 52 11

Introduction

The **52 11** is a pH electrode with screw cap S7, suitable for samples containing proteins. Large, semi-spheric membrane, very robust and fast responding.

Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it (5).

Verify that there are not air bubbles inside the membrane. They can be removed by shaking the electrode as a clinical thermometer. Remove the cup from the electrolyte filling port (3).

Calibration

 See pH-meter instruction manual.

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragms (7).
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

 Do not leave the electrode immersed in distilled water.

Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible (4).
- Maintain the membrane's surface clean (8). Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for HACH regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in HACH cleaning electrode solution or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Air bubble in the membrane.

Shake the electrode as a thermometer.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or substitute the electrode.

Storage

The connector should be protected from humidity by the protective cap ①.

Place the electrode in the storage protector ⑤ with KCl 3M + glycerine.

Close the refilling port with the stopper ③.

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not response appropriately after regeneration, must be replaced by a new one.

To avoid accidental breakage of the glass body electrodes, the use of the working protector is recommended.

 *If an electrode has too short lifespan, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

Measuring range	0... 12
Operating temp.	0... 100 °C
Reference element	encapsulated Ag/AgCl crystals
Diaphragms	3 ceramic
Electrolyte	KCl 3M + glycerine (LZW9502.99, LZW9512.99)
Body material	glass
Min. immersion depth	20 mm 

Applications

It is appropriate for samples with proteins and also for high temperatures up to 100 °C.

Limitations

Very viscous solutions or products with colloids able to clog the three diaphragms.

Warranty

HACH guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode **52 11** has passed the quality test and follows the next specifications:

- Asymmetry potential $< \pm 15$ mV.
- Sensitivity, pH 4...7 (at 25 °C) $> 98\%$.
- Response time, pH 4...7 < 20 s.

HACH COMPANY World Headquarters

P.O. Box 389, Loveland, CO 80539-0389 U.S.A.

Tel. (970) 669-3050

(800) 227-4224 (U.S.A. only)

Fax (970) 669-2932

orders@hach.com

www.hach.com

HACH LANGE GMBH

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D-40549 Düsseldorf, Germany

Tel. +49 (0) 2 11 52 88-320

Fax +49 (0) 2 11 52 88-210

info@de.hach.com

www.de.hach.com

HACH LANGE Sàrl

6, route de Compois

1222 Vérenaz

SWITZERLAND

Tel. +41 22 594 6400

Fax +41 22 594 6499